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Shaver blades unsteril

MD

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
Copyright by REMA Medizintechnik GmbH
All rights reserved
Technical changes reserved

Rev. V003

Date: 23.08.2022

ENGLISH

Caution

Please read all information contained in this insert. Incorrect handling and care as well as misuse can lead to premature wear or risks for the patient and user.

1 Scope

These operating instructions are applicable for REMA Shaver Blades non-sterile (product group 2AG4) risc class IIa:

REF 15-6XXXX – 15-79XXX.

2 Material

Tube: Stainless Steel 1.4301
Blade: Stainless Steel 1.4021
Top Part: Shaver Blades for
LINVATEC/Arthrex and Arthrex APS II - PMP
Dyonics EP1/Stryker SE-3/SE-4/SE-5 - Aluminium
Dyonics PS 3500 EP/Stryker Formula - PPSU
Aesculap/Arthronet - 1.4305

3 Intended use

The REMA Shaver Blades non-sterile are used for removal (cutting and milling) of delicate soft tissue or cartilage of small and large joints in arthroscopic surgery.

4 Indication

The following general indications belongs for the use of the arthroscopic treatment of the Shaver system:

- resection of meniscus tissue
- synovia resection
- resection of scar tissur
- aspiration of loose fragments and small meniscus particle
- cartilage debridement, removal of large cartilage fringes
- milling of exposed bones for induction of replacement cartilages
- notch plastic
- removal of osteophytes
- smoothing of drill holes entries

5 Contraindication

- local or generalized infections
- immune suppressive therapy
- coagulation disorders, embolism or thrombosis
- serious systemic disease that make a survival after the surgery unlikely
- if an arthroscopic treatment is not indicated in general

Product-specific contraindications are not existing at the moment.

6 Normal use

Only skilled medical personnel and personnel trained in this product are permitted to use Shaver Blades.

Please note before use:

Before initial use, new non-sterile Shaver Blades must be cleaned and also sterilized according to a validated procedure (DIN EN ISO 17665) (In this regard see Section "Reprocessing").

A visual test must be performed before each use (see section "visual inspection"). It must be ensured that Instrument-Insert is inserted firmly in the handle.

Application:

The type of appropriate attachments is selected depending on the indication and scope.

Blades bound for soft tissue should not be applied on the bone and vice versa.

Rema Shaver Blades have to be connected to shaver handpieces. They are used for the following shaver systems:

- REMA/ Linvatec/ /Arthrex
- Arthrex APS II/ REMA
- Atlantec/ Dyonics EP1
- Aesculap
- Stryker Formula

Please check before use if shaver blade is suitable for your system.

During use the shaver blades must be completely surrounded with irrigation fluid. Always ensure an adequate flow of irrigation fluid (at least 100 ml/min) to the joint.

Optimal Speed:

| Design: | rpm |
|---------|-------------|
| Blades | 800 – 2000 |
| Burrs | 2500 – 3500 |

Assembly and Operation:

- Position the shaver blade with the pen in the leadership of the handpiece.
- Shaver Blades are controlled by the foot pedal
- Please follow for introducing the REMA Shaver Blades in the Dyonics EP1 - handpiece the following instructions:

Assembly of Dyonics EP1 Shaver Blades

For insertion into the handpiece

1. Turn the knurled clockwise until it stops.
2. Blades in handpiece insert (pin into the small hole)
3. Turn the knurled counter-clockwise until stop to remove them from the handpiece
4. Turn the knurled right
5. Remove the blade from handpiece

ATTENTION:

- In case of **curved shaver blades** there must be determined the **position by a turning motion with a slight axial pressure** in which it is possible to insert or dismantle the blade effortless from the tube.

7 Safety notice

WARNINGS AND PRECAUTIONS!

- **Speed:** For the proper functioning of shaver blades meet the specified ranges of rotation. At too high speed, there is a risk that the removed material clogging the shaver. The speeds specified must be adhered to for safety and optimal use for. If anything is unclear, contact the manufacturer.

- **Damage:** If shaver blades operated in the dry state, this can result a damage (for example the blade).

The blades may in the rotating state don't come into contact with metal (eg, surgical instruments, arthroscope, etc.), otherwise it will cause damage.

The dismantling and the following assembling of the shaver blade (blade into tube) and the inserting of the shaver blade into the hand piece must be done carefully, in order to avoid damage to the shaver blades and / or injuries to the patient or surgical personnel. Shaver Blades may be damaged if excessive force is applied or by contact with metal!

Before you use the shaver blades they must be checked for damage. The Shaver Blade must be replaced by a new one (even during operation) if damage of any kind exists.

- **Cutting performance:** If a sufficient cutting power is not given, the Shaver Blades may no longer be used and must be replaced. See section "visual inspection".

8 Life cycle

The expected life cycle of the product depends of the proper use and the number of cycles.

Shaver Blades are by nature subject to increased wear depending on the type and time of use. Instruments-cutting geometry is dull with frequent use. The surgical manipulation cannot be performed satisfactorily in this case. Before each application of the product perform a visual inspection according section "visual inspection".

9 Storage / packaging / transportation / handling

Shaver blades must be stored in a dry location, protected from dust, moisture and sunlight at a moderate temperature in the range of 5°C to 40°C. They should be individually stored in a protective container with individual compartments or heat-sealed in film. The instruments must always be handled with the utmost care during transportation, cleaning, sterilization and storage. This applies in particular for fine tips and other sensitive areas. It is the operator's responsibility to ensure that the sterile condition is preserved after the sterilization process.

10 Reprocessing

Reprocessing cycles:

Due to the product design, the raw materials used and the intended purpose it is not possible to determine a precise limit with regard to the maximum possible number of reprocessing cycles. (See Section "Life cycle").

Preparation and transport:

The non-sterile shaver blades are to be taken from their packaging (for the first use). They must be placed in a designated for cleaning / sterilization container / device.

Disassembly of the instruments is necessary. For cleaning there must be separated (dismantled) the inner shaft from the outside shaft. The shaver blades should be sterilized pairwise after dismantling and after the reprocessing pairwise brought together.

REPROCESSING INSTRUCTIONS:

Preparation at the Point of Use:

Remove gross soiling by submerge the instrument into cold water (<40°C) immediately after use.

Don't use a fixating detergent or hot water (>40°C) as this can cause the fixation of residua which may influence the result of the reprocessing process.

Transportation:

Safe storage and transportation in a closed container to the reprocessing area to avoid any damage and contamination to the environment.

Preparation for Decontamination:

The devices must be reprocessed in an disassembled state.

AUTOMATED CLEANING PROCESS:

Pre-Cleaning:

Brushing:

Brush the disassembled instrument under cold water with a soft brush until all visible contamination is removed.

Rinsing:

Rinse shaver blades with the water jet pistol (static pressure above 4 bar) for a minimum time of 10 seconds.

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
Copyright by REMA Medizintechnik GmbH
 All rights reserved
 Technical changes reserved

Rev. V003

Date: 23.08.2022

ENGLISH

Washer/ Disinfec- Miele G 7735 CD (Miele)
 tor: Key Hole Surgery Rack E 450/1 (Miele)
 Detergent: Neodisher FA, Dr. Weigert
 Validation: SMP-Report-No.: 10109011407

Put the Shaver Blades in a disassembled state on a special key hole surgery rack and start the cycle. For couplings REMA, Linvatec, Arthrex, Arthrex APS II, Concept, Atlantec, Dyonics EP1 use for the outer tube and inner tube the REMA-Cleaning adapter.

| Step: | Time (min) | Process step | Reagents | Temp. (°C) |
|-------|------------|------------------------|--|------------|
| 1 | 1 | Pre-cleaning | Tap water | cold |
| 2 | | Drain | | |
| 3 | 3 | Pre-cleaning | Tap water | cold |
| 4 | | Drain | | |
| 5 | 5 | Cleaning | Tap water - 0,5 % alkaline detergent (Neodisher FA, Dr. Weigert) | 55 |
| 6 | | Drain | | |
| 7 | 3 | Neutralization | Tap water | cold |
| 8 | | Drain | | |
| 9 | 2 | Intermediate- rinse | Tap water | cold |

Automated Disinfection:

Automated Thermal Disinfection in washer/disinfector under consideration of national requirements in regards to A0-Value (see ISO 15883).

Drying:

Drying of outside of instrument through drying cycle of washer/disinfector. In-sufflate cavities of instruments by using sterile compressed air.

Visual inspection:

Perform visual inspection for cleanliness, perform an assembly and visual inspection according section "visual inspection". If necessary, repeat the reprocessing process until the instrument is optically clean.

Packaging (in case of packaging before sterilization):

Packaging has to comply with ISO 11607 and EN 868 standards for packaging for sterilised instruments.

STERILIZATION:

Sterilization of instruments by applying a fractionated pre-vacuum process (according. ISO 13060 / ISO17665) under consideration of the respective country requirements.

Parameters for the pre-vacuum cycle:

- 3 prevacuum phases with at least 60 milli bar pressure
- Heat up to a minimum sterilization temperature of 132°C; maximum temperature 137°C
- Minimum Holding time: 3 min (full cycle)
- Drying time: minimum 1min

Validation: SMP-Report-No.: 07213

Packaging (in case of packaging after sterilization):

Packaging has to comply with standards for packaging for sterilised instruments.

Storage:

Sterilised instruments have to be stored in a dry, clean and dust-free area at moderate temperatures from 5°C to 40°C.

MANUAL REPROCESSING:

Pre-Cleaning:

Immerse the instrument into cold tap water for at least 10 minutes.

Activate all movable parts several times up to the stop.

Brush the instrument with a soft brush until all visible impurities are removed.

Flush the Instrument with a water jet pistol for minimum 20 seconds.

Ultrasonics:

Immerse the instrument completely into an ultrasonic bath with enzymatic detergent (0,8%) and treat with ultrasonic (frequency 35 kHz) for at least 10 minutes at 45°C. Avoid sonic shadow.

Parameters:

- 45°C
- 10 min.
- 0,8 % cleaning solution
- 35 kHz

The instrument is taken out of the bath and flushed with a water jet pistol for minimum 20 seconds.

Then rinse the instrument with cold tap water again.

Ultrasonic bath: Bandelin Sonorex RK 1028 H

Detergent: Cidezyme/Enzol of company ASP (enzymatic)

Validation: SMP-Report-No.: 15812

Manual Chemical Disinfection:

Immerse the instrument for 10 minutes completely into a 4%-Mucocit-T solution (Company Merz Hygiene GmbH) at room temperature as specified by the manufacturer.

Validation: SMP-Report-No.: 26913

Drying:

Manual drying is carried out using sterile compressed air - in particular take care for drying cavities and channels.

Continuation of the process:

- Visual inspection
- PACKAGING
- STERILIZATION
- Storage

(as described in the above section "AUTOMATED CLEANING PROCESS".)

If the chemicals and machines described before are not available, the user is obliged to validate the process used.

11 Visual inspection

Before each use, the outer and inner shaft should be checked for pressure points or damage. Any damaged blades must not be used and replaced by new ones.

Especially the cutting edges should be inspected. They must be sharp and in good order. See section "safety notice – cutting performance".

12 Returns

Returns will be accepted only if they are marked as "hygienically safe" or "not decontaminated" and have been securely packaged for shipping.

Attention: Returned products must pass the complete reprocessing process before being returned for repair.

13 Liability / Warranty Claim / Repair and modification

REMA Medizintechnik GmbH as manufacturer and seller of these products assumes no liability for immediate or ensuing damages resulting from improper use or handling, especially due to failure to use as directed or improper processing or due to combination with products of other manufacturers.

Never attempt to perform repairs yourself. Service and repair work may only be performed by persons qualified and trained accordingly. Defective Shaver Blades may be sent to REMA Medizintechnik GmbH for repair. Unauthorized modifications and repairs are strictly prohibited and will entail invalidation of the manufacturer's warranty and liability.

14 Disposal

Shaver blades and the packaging material must be disposed of in accordance with the regulations and laws specific to the country in which they are used.

15 About these Instructions of Use

Throughout the period of use of medical devices, the Instructions for Use must be kept freely accessible for every user.

Explanation of the symbols:


| | |
|--|---|
| | Attention. Obey the documents accompanying the product regarding important safety relevant information as warnings and precautions. |
| | Medical Device |
| | Follow instructions for use |
| | Article number |
| | Batch code |
| | CE-Mark of the notified body mdc GmbH, Stuttgart |
| | Non-sterile |
| | Protect from moisture |
| | Keep away from sunlight |
| | Do not dispose of residual waste |
| | Manufacturer |

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