## ENGLISH

## INSTRUCTIONS FOR USE

## Retractors, self-holding

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Tel. +49(0)7461-92393-0  
Fax. +49(0)7461-92393-33  
www.dewimed.de

## IMPORTANT INFORMATION PLEASE READ BEFORE USE

Articles 18.xxx.xx; 40.xxx.xx; 54.xxx.xxx; 56.xxx.xx; 66.xxx.xxx, 68.xxx.xxx

In purchasing this instrument, you have acquired a superior-quality product, and information on its correct handling and use is detailed below. Please read and observe these user instructions carefully to minimise hazards for both patients and users in as far as possible.

Only skilled personnel who have been appropriately trained should be entrusted with use, disinfection, cleaning and sterilisation.

**Please note that our instruments are delivered in an unsterilised condition and should be cleaned, disinfected and sterilised prior to initial use.**

Self-retaining retractors and spreaders are instruments manufactured from metals such as steel, titanium and aluminium which are used in general surgery. Self-retaining retractors and spreaders are available in different models and versions and designed for short-term use (<30 days). For detailed information on available sizes and variants, please refer to our catalogue.

**Purpose / Indication:**

Self-retaining retractors and spreaders are used for spreading and holding back tissue, open wounds, bones, the thorax, sternum or other body orifices to enable a surgical intervention in the operative site.

Self-retaining retractors and spreaders are not intended for use in direct contact with the central nervous system or the correction of defects of the heart or central circulatory system!

Please note that the product LEMMON sternal approximator child is intended exclusively for use on children and patients of short stature. Selection of the appropriate spreader is incumbent upon the operator.

**Contraindications:**

- In general: Where a general operating risk exists in relation to patients.
- Self-retaining retractors and spreaders are not specifically intended for the monitoring, diagnosis, examination or correction of defects of the heart or central circulatory system in direct contact with these body parts.
- Self-retaining retractors and spreaders are not specifically intended for use in direct contact with the central nervous system.
- Self-retaining retractors and spreaders are not intended for the emission of energy in the form of ionised radiation.
- Self-retaining retractors and spreaders are not intended for triggering a biological effect or to be resorbed completely or to a significant degree.
- Self-retaining retractors and spreaders should not undergo any chemical change in the body.
- Self-retaining retractors and spreaders should not deliver any drugs. No other specific contraindications are known.

**User and safety instructions:**

Failure to observe these user and safety instructions may lead to injuries, malfunctions or other unexpected incidents.

All types of reusable self-retaining retractors and spreaders should be cleaned, disinfected and sterilised completely prior to initial use and any other use.

The self-retaining retractor or spreader must be inspected prior to every use for correct function, visible damage and wear (e.g. cracks or fractures). The function of self-retaining retractors and spreaders must be assured prior to every use.

The packaging is unsuitable for the high temperatures encountered during autoclaving and should be disposed of prior to initial sterilisation.

Do not overstress the instruments. Overstressing due to excessive application of force can cause breakage, bending and malfunctions in medical products and lead to injuries to the patient or user. Bent instruments should not be bent back into the initial position, as there is a risk of breakage.

Do not use a damaged or defective product. Withdraw any damaged product immediately, label and exclude from any further use.

**Reportable events:**

Please note: The laws in your country demand that, under certain circumstances, you report all serious incidents relating to use of our products directly after their occurrence to us, as the manufacturer, and/or our specialist dealer and to the responsible public authority.

**Use:**

Use of self-retaining retractors and spreaders is restricted exclusively to skilled medical personnel with special training.

Selection of a suitable self-retaining retractor and spreader is incumbent upon the experienced user.

**Following use:**

Self-retaining retractors and spreaders should be conditioned immediately after use.

**Handling of self-retaining retractors and spreaders:**

All surgical instruments should in all cases be handled with extreme care during transportation, cleaning, care, sterilisation and storage.

New instruments should undergo at least three machine cleaning cycles prior to initial sterilisation. This leads to the formation of a passivated layer on the surface that protects the instrument against discolouring and corrosion.

New instruments should be stored in indoor air without protective packaging in a closed cabinet / drawer. It is important to ensure that pertinent hygienic regulations are adhered to in each case. In the case of new instruments which are to be stored for a longer period of time, we recommend that they be removed from the sealed plastic bag and treated with a medical oil approved for sterilisation.

**General instructions for reconditioning:**

- Responsibility for professional treatment and preparation lies with the operator of the respective central sterile supply department and his or her assistants.
- The treatment and preparation of medical products which are to be used in a low bacteria or sterile condition should be realised employing suitably validated procedures and taking manufacturer specifications into consideration in a manner that guarantees the verifiable success of the procedure and ensures that the safety and health of patients, users or third parties are not endangered. Evidence of compliance should be compiled of treatment and preparation.
- Manufacturer specifications regarding the concentration, temperature, exposure, etc. of respective cleaning and disinfecting agents should be adhered to during treatment and preparation.
- Excessive concentrations of chemicals can lead to damage on instruments and contribute to rendering the laser or electrolytic labelling illegible.
- Treatment and preparation of instruments should be realised as soon as possible after use to avoid damage to the instruments.
- Professional treatment and preparation of surgical instruments begins on the operating table.
- New and repaired instruments should pass through the complete treatment and preparation cycle prior to use!
- It is imperative that instruments exhibiting visible damage be withdrawn immediately and/or professionally repaired by the manufacturer.
- Damaged instruments can jeopardise the success of an operation!

**Process for treatment and preparation of self-retaining retractors and spreaders:****Preparation at location of use:**

Directly after use, coarse dirt should be removed from instruments by immersion in cold tap water (<40 °C). Do not use any fixing detergents or hot water (>40 °C), as this will lead to fixation of residue which may negatively influence the success of cleaning.

**Transportation:**

Instruments should be kept in a closed container for safe transportation. This should prevent any damage to the instruments

**Preliminary cleaning:**

- The self-retaining retractors and spreaders are disassembled in as far as possible and the components immersed in cold tap water for 10 minutes, ensuring that all component parts of the self-retaining retractors and spreaders are covered with water.
- Using a soft brush, thoroughly clean all component parts of the self-retaining retractors and spreaders individually under flowing cold tap water until all visible impurities are removed.
- Repeat this process until all visible impurities have been removed.
- Subsequent to this, thoroughly rinse all positions difficult to access such as hinges, contact surfaces, cavities, bored holes and threads with the water pistol (water jet pistol with a static water pressure of >2 bar) for a minimum of 20 seconds.

In the case of severe dirt or caking, preliminary cleaning is recommended in a US bath with the following parameters:

- Place the instruments in the metal sieve so that they do not touch each other.
- Immerse completely in a 0.8% Cidezyme solution.
- Clean at room temperature and 35 kHz for >10 min.
- Following preliminary ultrasound cleaning, the instruments are removed and cavities, bored holes and threads rinsed with cold tap water from a water pistol for a minimum of 20 seconds in pulse mode.

**Manual cleaning:**

Manual cleaning can only be recommended if machine conditioning in a cleaning and disinfecting unit is impossible.

- Place the instruments in the metal sieve so that they do not touch each other.
- Immerse completely in a 0.8% Cidezyme/Enzol or Mucadont Zymaktiv cleaning solution.
- Clean at 45 °C and 35 kHz for >10 minutes.
- Following ultrasound cleaning, the instruments are removed and cavities, bored holes and threads rinsed with cold tap water from a water pistol for a minimum of 20 seconds in pulse mode.
- The self-retaining retractors and spreaders are finally rinsed with cold, demineralised water.

**Manual disinfection:**

Manual disinfection can only be recommended if machine conditioning in a cleaning and disinfecting unit is impossible.

The instruments are immersed completely in a 4% Mucocit-T solution and disinfected at room temperature for between 5-10 minutes.

After removal from the disinfection bath, lay the instruments in demineralised water and rinse thoroughly. Dry after rinsing and realise sterilisation with damp heat.

**Machine cleaning, thermal disinfection and drying:**

Machine processes were validated with a Miele cleaning and disinfecting unit, model 7836CD. Place the self-retaining retractors and spreaders into a mesh tray on the infeed carriage. Start the cleaning process with the following minimum settings:

- 2 min. preliminary cleaning with cold tap water (16 °C ±2 °C)
- Drain
- 5 min. cleaning at 55 °C. Dispense 0.5% MedClean forte with tap water
- Drain
- 3 min. neutralisation with cold demineralised water (20 °C ±2 °C)
- Drain
- 3 min. rinsing with cold demineralised water (20 °C ±2 °C)
- Drain

**Settings for thermal disinfection:**

Realise machine thermal disinfection, taking national requirements regarding the A0 value into consideration (see ISO 15883).

- 2 min. rinsing with warm demineralised water (>40 °C)
  - Heat to disinfecting temperature 293 °C\*
  - Holding time at >90 °C\* for ≥10 min.
- \*Disinfecting temperatures relate to the upper and lower switching points of the cleaning and disinfecting unit thermostat

**Drying:**

- The cleaning and disinfecting unit program should assure a drying duration of at least 20 minutes at max. 93 °C.
- Following expiry of the drying duration, the instruments are immediately taken out of the cleaning and disinfecting unit.

**Manual drying:**

If necessary, additional manual drying can be realised with the aid of a lint-free cloth. Dry cavities with sterile compressed air.

**Care, functional testing and packaging:**

A careful visual inspection of the self-retaining retractors and spreaders for cleanliness is conducted following the cleaning and disinfecting process. If residual impurities are identifiable, the entire cleaning and disinfecting process should be repeated until residue-free cleaning can be confirmed.

If cleaning is not possible, the instrument should be withdrawn and excluded from further use. The instrument should be disposed of in this case.

The instruments and their component parts should then be inspected for possible damage such as cracks and concealed, loose or missing parts. Reassemble all component parts subsequent to this and check the function of the instrument.

Treat threads and hinges with an approved medical maintenance oil. The oil used (e.g. paraffin conforming to Ph. Eur.) should not affect subsequent sterilisation.

Withdraw defective instruments, label them and exclude from any further use. Take note of any risk of injury if appropriate.

**Packaging**

Only use standardised and approved packaging materials and systems (EN 868 Part 2-10, ISO 11607 Part 1 + 2, DIN 58953).

The self-retaining retractors and spreaders should be packaged in a manner that excludes the possibility of damage to the sterile barrier.

Instruments should not be sterilised in the protective and transportation packaging sent with the delivery

**Service lifespan:**

The service lifespan of surgical instruments is only insignificantly affected by the number of conditioning cycles realised if they are conducted in accordance with the validated procedures described here. It depends to a greater degree on considerate and careful handling of the instruments during all usage, conditioning, transportation and storage phases. The end of the service life is reached if evidence of wear or defects which limit the function of the product are detected during the prescribed visual inspection and functional testing. It is imperative that the instruments be marked in this case, excluded from any further use and replaced by functional instruments. Furthermore, the end of the usage cycle is reached if clear identification of the instruments is no longer possible due to the absence of the marking. It is imperative that instruments exhibiting visible damage be withdrawn immediately and/or professionally repaired by the manufacturer.

Please note: If disposable instruments are mixed with reusable surgical instruments during treatment, the entire contents of the cleaning and disinfecting unit may suffer corrosion, or the entire operation tray may rust in the steam sterilizer.

In all cases, we recommend machine cleaning and thermal disinfection of surgical instruments prior to sterilisation with damp heat..

**Sterilisation with damp heat:**

Sterilisation of products is preferably realised using the fractionated pre-vacuum method (pursuant to ISO 13060 / ISO 17665 and EN ISO 285) while taking respective national requirements into consideration.

- 3 pre-vacuum phases with at least 60 mbar pressure
- Heating to a sterilisation temperature of typically 134 °C
- Holding time: typically 5 minutes
- Drying time: minimum 10 minutes

The above specifications are recommendations of RKI/KRINKO for sterilisation with damp heat. These specifications were validated for our self-retaining retractors and spreaders with the SMP No. 23616 study at reduced settings (132 °C, 4 min. holding time).

Other regional and national regulations and directives may apply.

With regard to the treatment and preparation of medical products, a separate conditioning process should be employed in the case of patients suffering or suspected of suffering from Creutzfeldt-Jacob disease (CID) or its variant (vCID), or the instrument should be disposed of.

**Storage:**

The storage conditions of the packaging manufacturer for maintaining an effective sterile barrier apply when it comes to the storage of sterilised instruments. The instruments themselves do not require any special storage conditions.

**Repairs:**

Do not carry out unilateral repairs. Servicing and repairs should only be realised by suitably qualified personnel. Please contact the manufacturer or your medical technology department should you have any questions in this respect.

Defective products should undergo the entire reconditioning process prior to being returned for repair. Moreover, a Hygiene Clearance Certificate should be included in each return delivery.

**Disposal:**

No special measures are necessary with regard to disposal. Observe respectively applicable national or international laws and regulations during disposal.

**Treatment validation, studies, information:**

The following materials and machinery were used to validate individual treatment steps. Information in this respect can be obtained from the manufacturer.

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Machine cleaning: SMP No. 16016

Sterilisation: SMP No. 23616

It is the responsibility of the user to validate its process appropriately if the previously described chemicals and machinery are not available.




It is the obligation of the user to ensure that the reconditioning process and the necessary operating media, materials and personnel are suitable for achieving the required results.

The technological state of the art and national laws demand that these processes and the operating media used are kept in a validated and maintained condition.

**Explanation of symbols**

	Article number		Lot-number		Follow instruction for use
--	----------------	--	------------	--	----------------------------

	Medical device		Manufacturer		Non-sterile
---	----------------	---	--------------	---	-------------

	Quantity in packaging		Manufacturing date		Attention
---	-----------------------	---	--------------------	---	-----------

	Prescriptive medical device (USA)		CE-Mark of the Notified Body mdc medical device certification GmbH, Stuttgart, Germany
---	-----------------------------------	---	--

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Baño ultrasónico: Merz Hygiene GmbH: Mucadont Zymaktiv
Desinfectante: Merz Hygiene GmbH: Mucocit-T
Detergente: (Chemische Fabrik Dr. Weigert GmbH & Co. KG, REF 405033)