- During exposure time fill a 20-ml syringe with disinfectant and flush the hinge and all hard-to-reach and accessible areas of the device for at least 5 times. Pay attention to the critical and hard-to-reach areas where it is not possible to assess the cleaning efficacy for the devices by means of visual inspection.
- After 5 minutes exposure time start ultrasonic for 5 minutes.
- Turn off ultrasonic.
- Remove the device from the disinfectant.
- Rinse the device under running tap water for at least 1 minute.
- Repeat the rinsing step with deionized water.
- Wipe the device with single use, fluff free wipe or towel or dry by an air gun with compressed air.

Cleaning and Disinfection - Automated

- If visible soil remains in gaps, holes, joining surfaces, threads and other hard to reach areas: Start with Precleaning (see pre-cleaning above) and continue with manual cleaning.
 - Equipment:** Washer disinfecter in accordance with DIN EN ISO 15883-1+2 with thermal program (temperature 90-93oC), mildly alkaline cleaner detergent (e.g. Sekumatic® MultiClean 200 ml), fluff-free wipe or towel or air gun connected to compressed air of medical grade (in accordance with European Pharmacopeia).
- Bring all parts of the devices into the washer disinfecter and adapt them with their cannulated lumens to a load carrier module for minimally invasive surgery (MIS).
 - After closing the door start the thermal program, the program parameters are shown in the table below.
 - At the end of the program remove the device.
 - Check if the device is completely dry, if necessary dry with single use, fluff free wipe or towel or dry by an air gun with compressed air.

Program Step	Water	Dosage	Time	Temperature
Pre-Rinse	Cold		5 Min.	
Dosage cleaner		According to manufacturer's instructions		According to manufacturer's instructions
Cleaning	Deionized water		10 Min.	55°C
Rinse	Deionized water		2 Min.	
Disinfection	Deionized water			Ao-Wert >3000 ¹ e.g.: 5 Min. 90 °C
Drying			15 Min.	Up to 120°C

¹ Public authorities may adopt other implementing provisions (parameters for disinfection performance) within their field of responsibility

Maintenance, Inspection and Testing

- Equipment:** Silicone-free steam sterilizable lubricant allowed for use on medical products according to the manufacturer (e.g. Sterilit®, spray of B. BRAUN Aesculap or other fabricates)
- Inspect all components of the device under normal lightning for the removal of soil. Repeat manual cleaning (including pre-cleaning) and automated cleaning and disinfection or manual disinfection if components are not visible clean.
 - Check the device for completeness, damage and functionality.
 - Check the device with long slender features for distortion.
 - Do not use damaged instruments.
 - Lubricate the moveable instrument parts using silicone-free steam sterilization lubricant.

Packaging

Packaging in sterile barrier system in accordance with DIN EN ISO 11607 or DIN 58953. (Double packed sterilization paper wraps were used for validation

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- Pre-vacuum process, 134oC, and sterilization time at least 3 minutes (longer holding times and temperatures (<141oC) are possible).
 - Sterilization of the device in the disassembled status
- Place the package device in the sterilizer chamber.
 - Start the program.
 - At the end of the program remove the device and let it cool down.
 - Check that the seal of the package is closed, and the package is dry.

Storage

Storage and shelf life as established by the processor.

- Instruments should be stored in dry, clean and humidity-free environment without direct solar irradiations.
- Do not store them in places where chemicals may emit corrosive fumes or where fluctuating temperature of humidity may cause vapors to condense on the instruments.

Additional Information

Only validated processes may be used for the processing of medical devices

Symbols used on Labeling

	Manufacturer		Catalogue Number		Lot Number
	Warning		Non Steril		Consult Instructions for Use
	Keep dry		Protect from sunlight		
	CE marking and identification number of notified body mdc medical device certification GmbH, Stuttgart, Germany				

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