
ESPAÑOL		
INSTRUCCIONES DE USO		Electrodos artroscópicos
	Dewimed Medizintechnik GmbH Unter Hasslen 14 78532 Tuttlingen GERMANY	Tel. +49(0)7461-92393-0 Fax. +49(0)7461-92393-33 www.dewimed.de

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Sterilization

Only products that have been cleaned and disinfected are allowed to be sterilized. Only apply steam sterilization in an autoclave (fractioned pre-vacuum with sufficient product drying) for this product.

Adjust sterilization parameters:

- Minimum 134°C and maximum 137°C in saturated steam.
- Holding time at least 5 Min. until max. 20 Min.
- Drying in vacuum for at least 10 Min.

Example: Sterilizer Class B, manufacturer: Tuttnauer.

Respect data of the sterilizer manufacturer regarding the load, handling, and drying times.

Exclusion:

Do not apply hot air, EO-gas, Radiation, or Plasma for sterilization, or any other sterilization method for this product.

IMPORTANT:

Before usage, let the product cool to room temperature. It is the sole responsibility of the user to maintain the sterile condition of the product after the sterilization process. In case of the aforementioned chemicals and machines for cleaning, disinfection, or sterilization are not available, it is the responsibility of the user to validate the procedure actually applied. Also, if a sterilization method other than that described above is applied, this deviating procedure has to be validated by the user accordingly.

Limitation of Reprocessing

Product life depends on wear and tear, handling, application time, damage and frequency and type of reprocessing. Due to the materials and the intended use, this product may only undergo a maximum of up to 5 reprocessing cycles (see section "Visual and functional inspection"). Therefore, a visual and functional check must be carried out before each reuse. Only an undamaged product may be reused.

Visual inspection and Function Test

Before usage, the complete product, especially insulation and product tip have to be inspected for pressure points and damages. Products exhibiting damage or pressure points are not allowed to be used and have to be replaced with new ones. During and after application, tissue may adhere to the product, or sooting may be present on the distal end of the active electrode. Such adhesions or sooting do not represent a reason for complaint and the product has to be exchanged for a new one. Due to longer application time, mechanical forces, plasma seam, or similar, the product may exhibit deformation or abrasion of the insulation material. Also, such aspects do not represent a reason for complaint and the product has to be replaced with a new one. Blockage of the suction channel (if applicable), does not represent a reason for complaint. A product, that exhibits a blocked suction channel (if applicable), has to be replaced with a new one. Before usage, an electrical continuity test has to be done. In case the product does not pass the electrical continuity test, it is not allowed to use the product any longer and has to be replaced with a new one. If damaged, do not use these products. Never lay this product on the patient or in their direct vicinity.

Exclusion of Repair and Modification

Defective products are not allowed to be repaired. Unauthorized modification or repair is strictly prohibited and leads to an immediate loss of the manufacturer's warranty. It is not allowed to bend products or manipulate the product in any other manner. Especially products with a hook as an active part are never allowed to be bent. This could lead to severe injuries to the patient, the user, or the third party. The product is supplied non-sterile and must be cleaned, disinfected and sterilized before use according to the validated procedure (DIN EN ISO 17665) specified in these instructions for use. Before reuse, a visual inspection and function test must be carried out.

Only products that have passed the visual inspection and function test may be used again. In particular, undamaged insulation and an undamaged electrode tip are decisive for whether a product may be used again or not.

Packaging, Storage, Transport

The product must be stored in a clean and dry environment, they must be protected from direct sunlight. Always handle the product with the utmost care when transporting, cleaning, disinfecting, maintaining, sterilizing, and storing it. This is especially true for fine tips and other sensitive areas (e.g., insulation). Do not store or transport the product together with sharp-edged or pointed objects. Store only in protective containers with individual compartments or individually shrink-wrapped. The maintenance of the sterile state after the sterilization process must be ensured by the operator.

Safety, installation and application

Only to be used by skilled medical professionals who have been introduced to the product. In addition to the acknowledged benefits of HF technology, the application includes several risks that have to be attended to:

Improper use and failure to follow the instructions for use may result in unintentional burns to the patient, user, or third parties. Continued further education of the surgical personnel is recommended. It is the operator's responsibility to ensure the training of users/surgeons in arthroscopic surgical techniques. The HF generator, which is used in combination with the product, has to be checked for damages. Follow the instructions outlined in the user manual of the manufacturer.

Preparation:

Prepare the HF generator according to the User Manual and place it for application. Proceed with the generator according to the User Manual. Bring the sterilized product into the sterile field and confirm the integrity of the product's insulation by visual inspection. Confirm the integrity of the insulation of the product by visual inspection. At the beginning of the application, make sure that the cable of the handle is connected properly to the HF generator and that the correct settings are selected and displayed. It is the responsibility of the user to select the lowest possible power setting on the HF generator to achieve the desired effect during the respective operation.

Note: when using products with suction, connect the suction adapter to a tube that leads to a suitable suction system (e.g., suction canister, wall suction system, etc.). The minimum vacuum target value is 300 mmHg.

disinfectants if applicable. Thoroughly remove the surface stain with a soft brush or synthetic fleece, as otherwise particles or dried secretions may adhere. Ensure that hard-to-reach areas are well cleaned and rinsed several times. This may make subsequent cleaning and sterilisation difficult or impossible. Cavities and lumen have to be rinsed intensively using at least 3x20 ml cold tap water (<40°C) with the aid of a rinsing adapter (e.g., from the company Medisafe), or with a syringe or with a water jet pistol (>30 Sec.). This pre-cleaning step always has to be done prior to the manual cleaning or cleaning with the cleaning and disinfection device (CDD).

Manual Cleaning and Disinfection

Prepare an immersion bath with a suitable fluid cleaning agent. Use a cleaning agent that is compatible with the disinfectant and suitable for immersion baths. Follow the instructions of the manufacturer of the cleaning agent and disinfectant. Only use agents suitable for medical devices made from metal and plastics with a pH value between 5.5 and 12.3.

- **Recommendation:** Cleaning agent gigazyme® (Schülke & Mayr) and disinfectant Korsorex Plus. Do not use high alkaline cleaning agents. These will impair the lifetime of the product.
- Prepare immersion bath with cleaning agent according to the specific cleaning agent IFU.
- Prepare a separate immersion bath with disinfectant according to the specific disinfectant IFU.
- Immerse the product completely in the ultrasonic bath with a cleaning agent (e.g., 0.5% gigazyme®)
- Clean product in an ultrasonic bath using a sonication time of 5 Min. and a frequency of 45 kHz.
- Follow all instructions outlined in the IFU of the cleaning agent, disinfectant, and ultrasonic bath.
- Ensure the product will not touch other products or parts in the ultrasonic bath.
- Ensure sonic shadows in the ultrasonic bath are avoided.
- Then, clean the product with a soft brush under cold running town water (<40°C).
- Intensively rinse cavities and lumen with a water jet pistol (>30 Sec.) or similar for at least 1 Min.
- Afterwards, rinse the product thoroughly for at least 1 min. under running tap water (>40°C) to remove any residues of the cleaning agent.
- Inspect the product visually on the remaining stain.
- In case stain is still present, repeat the aforementioned cleaning steps as long as it needs until no soil is present.
- Afterwards: immerse product completely in a disinfectant bath including e.g., Korsorex Plus 3%, for at least 15 Min.
- Follow the manufacturer's data for residence time.
- Ensure the disinfectant will contact all areas on the product.
- Rinse cavities and lumen several times, which means at least 3x with 20 ml each of disinfectant bath fluid.
- Afterwards: Rinse the product thoroughly for at least 1 Min. with demineralized cold water, to remove all disinfectant residues.
- Additionally: Rinse all narrow and areas difficult to access, all cavities and lumen with a syringe several times (at least 3x) using each time 20 ml cold demineralized water.
- Dry product with a lint-free wipe and sterile compressed air.
- Dry cavities, lumen, and channels with sterile compressed air.

Automated Cleaning and Disinfection

Only use cleaning and disinfection devices (CDD) with proven efficiency according to EN ISO 15883. Follow the data of the manufacturer of the cleaning and disinfectant machine. Only use agents suitable for medical devices made from metal and plastics with a pH value between 5.5 and 12.3.

Recommendation: neodisher® mediclean forte (Dr. Weigert GmbH & Co. KG).

Apply program for thermal disinfection. Follow instructions and data regarding program courses and machines. Product has to be stored safely and protected against mechanical damages during automated cleaning and disinfection. Do not clean together with sharp-edged or pointed objects. Deposit product in a suitable rinsing basket. Follow data for loading of the cleaning and disinfection device (CDD). Use rinsing adapters for products with the lumen and connect them according to the instructions in the User Manual of the cleaning and disinfection device (CDD).

Cleaning Program

Start the program course with the following parameters:

- 1 Min. pre-rinsing with cold water
- Emptying
- 3 Min. pre-rinsing with cold water
- Emptying
- 5 Min. cleaning at 55°C with 0.5% alkaline cleaning agent
- Emptying
- 3 Min. neutralization with warm tap water (>40°C) and neutralizer (0.1% Neodisher® 2)
- Emptying
- 2 Min. interim rinsing with warm demineralized water (>40°C)
- Emptying

Disinfection Program

Automated thermal disinfection considering national requirements regarding AO value (see EN ISO 15883, AO value >3000).

- 5 Min. cleaning at 92°C +/-2°C

Drying

- 30 Min. at 90°C
- Remove rinsing adapter

At the end of the program course, remove the product and inspect it on the remaining stain. In the case of residues, repeat the automated cleaning and disinfection step as long as it takes until the stain is no longer present.

Dry cavities and not sufficiently dried areas with sterile compressed air <2 bar. Immediately after removal of the product and immediately after additional drying in a clean place, put the product in a single-use sterilization packing (double packing) from paper or foil or put the product in a sterilization container. Respect requirements for sterilization packaging according to EN ISO 11607 and EN 868.

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