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**HF-Kabel und Adapter**

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**Scope**

Maximum rated voltage of accessory (Umax):

Description	Umax
Monopolar Kabel/Adapter	4,3 kVp
Bipolar Kabel / Adapter	0,8 kVp
Neutraelektrodenkabel	4,3 kVp

**See also label or catalogue.**

In any combination with another electro-surgical accessory, the maximum rated voltage of the combination corresponds to the lowest rated voltage of the accessories used. (See also section "Intended Use"). In case of uncertainties, contact the manufacturer. Read the entire IFU of this product and any used accessory as well as HF generator and HF neutral electrode (monopolar application). All requirements, safety notices, and warnings included in the respective IFUs have to be followed strictly.

**Intended Use**

**The product must be connected to the designated output of the electro-surgical generator and serves as a contact element between the HF generator and monopolar and bipolar electrodes. Use only by trained medical personnel.**

The application of HF current may cause damage to pacemakers or in vivo cardiac defibrillators, therefore affected patients must consult a cardiologist before the procedure.

**Connection and Activation:**

Before using monopolar products, ensure that a designated HF neutral electrode is correctly applied to the patient and properly connected to the designated HF generator. Connect product only to switched-off HF generator or in standby mode (where applicable).

**Compatibility:**

Before use, ensure that the product is compatible with the HF accessories and HF generator intended for the procedure.

In case of uncertainty, contact the manufacturer of this product or the manufacturer of the accessories, HF generator, etc. used.

Use only with compatible products and connectors.

Never connect to different or unknown power sources!

Using the adapter, activation can only be performed via the foot switch of the generator, using the bipolar HF cable.

The instructions, safety and warning notes of the operating instructions of the accessories used and the HF generator used must be followed.

It is recommended to use a smoke evacuation system.

**IMPORTANT:**

The product must be handled with the utmost care.

This applies not only for the duration of the procedure, but also for the duration of storage, preparation, storage and transport, as well as during the process of connecting to the HF accessories and to the HF generator.

This applies in particular to the fine components and other sensitive areas, such as the insulation.

Improper use will result in immediate loss of warranty.

No liability is accepted for any damage that may occur.

**Safety Notices – WARNING**

The maximum rated voltage of the product can be found in these operating instructions, on the label or in the current product catalog.

In case of uncertainty, contact the manufacturer.

The product must be cleaned and disinfected as well as sterilized before first use and after each use according to a validated procedure (DIN EN ISO 17665) (see section "(Re)processing:: Cleaning, disinfection and sterilization").

A visual and functional check must be performed before each use (see section "Visual and functional check").

If damage, deformation or similar is found on the product, the product must not be used.

It must be replaced by a new product.

At least one (1) cleaned, disinfected and sterilized replacement product must be kept on hand. It is the user's responsibility to determine the appropriate product size and type based on professional judgment and the patient's specific indication, preferred surgical technique and history, etc.

Before use, ensure that the product is firmly and correctly connected to the HF generator and that the HF accessories used are firmly and correctly inserted.

Care must be taken to avoid damage to the product and/or injury to the patient, surgical staff or third parties.

The product can be damaged by excessive force.

The product must therefore be observed throughout use.

Moisture in the connectors can cause voltage flashovers which damage the connector or its insulation.

Avoid penetration of liquids into the plug during use.

**Exclusion:**

Only grasp the product by the plug part when plugging in and unplugging it.

Do not activate the HF generator while the product is connected.

Do not bend the product.

Do not use the product coiled or rolled up.

Do not activate the product while it is in contact with metallic objects and/or optics.

During an electro-surgical procedure, the patient must not come into contact with grounded metal objects such as a surgical table frame, instrument tray or similar.

Ensure that there are no flammable substances in the immediate vicinity during the entire application, otherwise there is a risk of explosion.

**(See also the section entitled "General safety Instructions for using HF technology").**

After the electro-surgical current has been switched off, the product tip may still be hot so that it can cause unintentional burns.

The product must not be activated continuously for a long period of time.

It is the user's responsibility to select the lowest possible power setting on the HF generator in order to achieve the desired effect during the respective procedure.

During the surgical procedure, mechanical forces may cause deformation or wear of the product.

The application of RF current may cause damage to pacemakers and in vivo cardiac defibrillators, therefore affected patients must consult a cardiologist before the procedure.

**General Safety Notices for HF-Technology (Excerpt)**

In addition to the acknowledged benefits of HF-Technology, the application includes several risks that have to be attended to:

Improper use and non-observance of the IFU can lead to unintentional burns on the patient, user, or third party. Continued further education of the surgical personnel is recommended.

**a) Environment**

Pay attention that no flammable substances (anesthetics, oxidizing gases, endogenous gases, etc.) are present in the immediate vicinity as otherwise there is a danger of explosion.

Only use non-flammable cleaning or disinfection agents.

All oxygen connections must be tight and leak-proof during the procedure.

**b) Patient Positioning and Preparation**

Ensure proper patient positioning, i.e., only use insulating surgical drapes that are dry, absorbent, and liquid-tight. Isolate conductive surfaces and contact points towards the patient. Dry tissues of pulp are necessary to use for skin and breast folds as well as between the extremities.

Prior to application, remove any liquids that potentially accumulated in body cavities.

Only use non-flammable disinfectants. Do not use, for example, alcohol-based tinctures.

Only use non-conductive rinsing fluid, if medically possible. Attend the requirements on irrigation fluid for monopolar and bipolar products. Prior to application, remove any body jewelry from the patient. Putting a band-aid over the body jewelry is not sufficient!

**c) Connections**

Prior to application, make sure the product is connected properly to the HF-generator and make sure the correct power and performance setting is adjusted and displayed. It is the responsibility of the user to apply low-performance settings on the HF-Generator to achieve the envisaged effect for the respective intervention.

**d) HF-Neutral Electrode for Monopolar Application**

In the case of monopolar application, select an HF neutral electrode suitable for the patient, apply it correctly and connect it properly with the respective HF generator.

Follow all instructions for proper application of the HF neutral electrode, incl. patient protection and patient monitoring, monitoring of the HF neutral electrode, and all further provisions, safety notices, and warnings included in the IFU of the HF-neutral electrode.

**e) Patient Reactions**

All electro-surgical devices potentially can cause muscle stimulation during the application. The design of this product minimizes the risk of this undesirable effect. Nevertheless, muscle stimulation can lead to an unexpected movement of the patient in the surgical field.

**f) Handling HF-Accessory**

Make sure the accessory is compatible. Do not touch the instrument tip during the complete application. As long as the product is not applied, place it on a dry, clean and non-conductive and well-visible surface that is not in contact with the patient. Never store product on the patient. Unintended activation of the product can lead to unintentional burns or other injuries to the patient, the user, or third party.

Never wrap cables around the patient and never lay cables over the patient. Apply only short activation times. Keep on longer breaks between activation phases. Only adjust low power settings.

**g) Completeness of the System**

At the end of the surgery, confirm the completeness of the system.

**(Re) Processing, Cleaning, Disinfection and Sterilization**

**Maximum number of reprocessing cycles.**

**Due to the design, materials used, intended use as well as wear and tear, a maximum limit of performable reprocessing cycles cannot be determined.**

When applied according to the Intended Use, the product underlies natural wear and tear, considering the manner and duration of the application as well as the manner and frequency of reprocessing.

Therefore, a visual inspection and function test has to be done prior to each usage. (Refer to section "visual inspection" and "function test"). Visual inspection and function test, especially the condition of the complete cable, insulation, port, plug component, and pins are decisive for whether the product is allowed to be applied again.

**Time requirements for cleaning and disinfection.**

Preparation for cleaning, pre-cleaning, and automated cleaning and disinfection has to be done immediately after the application, however not later than 1 hour after the application. Avoid idle time.

**After the Application**

Clean and disinfect the product immediately after the application. However not later than 1 h after application. After the application of the product, deposit the product carefully (protection of the life of the product). After application, separate the contaminated product and deposit it in a suitable container (deposit means "do not drop"). Immediately remove the gross stain. Immediately mark damaged or defective products.

Accessory, that does not fit the sieves of the cleaning and disinfection device (CDD), shall be deposited separately in suitable containers. Close firmly all disposal-/transport containers, to avoid drying of stain. Organize transport of contaminated products the way that contamination of the transport ways and environment also is avoided (closed transport).

Unused reusable products have to be reprocessed as well.

Take care, all and any transport containers are cleaned and disinfected after transport as well.

**Validation of (Re)Processing**

The following validated processing procedure is recommended. Equivalent deviant processes are possible. Then, it is the sole responsibility of the user to safeguard the suitability of the actually applied procedure by suitable means (e.g., validation, routine monitoring, verification of material compatibility, etc.). Automated cleaning and disinfection always are preferable. The following procedure was validated according to EN ISO 17665.

Additionally, applicable processing requirements specific to the respective clinical place (operator) as well as national or country-specific regulations have to be followed as well.

Never use sharp objects for cleaning. Disinfectants always have to be rinsed carefully.

**Preparation for Cleaning**

New product: Remove the product from its packaging.

Used product: Any attached HF accessory (e.g., HF-Adapter, HF-Electrode, etc.) has to be disconnected. Place the product in a container intended for pre-cleaning.

**Pre-Cleaning**

The product has to be cleaned immediately after usage (see section above).

Prepare a cleaning bath using an 8 ml cleaning agent (e.g., CIDENZYME®) per 1 l of water.

Attend requirements set forth by the manufacturer of the cleaning agent.

**Requirements for cleaning agent:**

Intended for medical devices made from metal and plastics.

pH value: 5.5 to 12.3

Do not use high alkaline cleaning agents. These will impair the lifetime of the product.

Attend the reaction time set forth by the manufacturer of the cleaning agent. Take care the product will not touch other parts in the cleaning bath. Thoroughly remove the surface stain with a soft brush or synthetic fleece, as otherwise particles or dried secretions may adhere. This could make subsequent cleaning and sterilization difficult or impossible. Do not use sharp, sharp-edged, or scrubbing devices for cleaning. Ensure that areas difficult to reach are cleaned thoroughly and rinsed several times. Rinse the product thoroughly for at least 1 minute under running tap water, to remove residues of the cleaning agent.

Check the product visually for remaining stains. In case the stain is still visible, repeat the process of pre-cleaning as often and long as it needs until all stain is removed and stain is no longer visible. This pre-cleaning step always has to be done prior to the manual cleaning or cleaning with the cleaning and disinfection device (CDD).

**Manual Cleaning and Disinfection**

Manual cleaning and disinfection are not recommended.

**Automated Cleaning and Disinfection**

Only use cleaning and disinfection devices (CDD) with proven efficiency according to EN ISO 15883. Follow the data of the manufacturer of the cleaning and disinfectant machine.

**Requirements for cleaning agent:**

Intended for medical devices made from metal and plastics.

pH value: 5.5 to 12.3

Do not use high alkaline cleaning agents.

These will impair the lifetime of the product.

Recommendation: neodisher® medclean forte (Dr. Weigert GmbH & Co. KG).

**Requirements for CDD:**

Apply program for thermal disinfection.

Follow instructions and data regarding program courses and machines.

Product has to be stored safely and protected against mechanical damages during automated cleaning and disinfection.

Do not clean together with sharp-edged or pointed objects. Deposit product in a suitable rinsing basket. Take care cable is not kinked or squeezed. Follow data for loading of the cleaning and disinfection device (CDD).

**Cleaning and Disinfection Program:**

Start the program course with the following features:

Thermal disinfection

5 to 10 minutes at 90°C to 93°C (tolerance according to EN ISO 15883-1 or -2) A0 ≥ 3000

Final rinsing with deionized water

Product drying.

At the end of the program course, remove the product and inspect it on the remaining stain.

In the case of residues, repeat the automated cleaning and disinfection step as long as it takes until the stain is no longer present. Dry cavities and not sufficiently dried areas with sterile compressed air <2 bar. Immediately after removal of the product and immediately after additional drying in a clean place, put the product in a single-use sterilization packing (double packing) from paper or foil or put the product in a sterilization container.

Respect requirements for sterilization packaging according to EN ISO 11607 and EN 858.

**Sterilization**

Only products that have been cleaned and disinfected are allowed to be sterilized.

Only apply steam sterilization in an autoclave (fractioned pre-vacuum with sufficient product drying) for this product.

**Adjust sterilization parameters:**

Minimum 134°C and maximum 137°C in saturated steam.

Holding time at least 5 Min. until max. 20 Min.

Drying in vacuum for at least 10 Min.

Sterilizer (Class B) according to valid national standards and regulations (e.g., EN 13060 or EN 285).

**Example:** Sterilizer Class B, manufacturer: Tuttnauer.

Respect data of the sterilizer's manufacturer regarding the load, handling, and drying times.

**Exclusion:**

Do not apply hot air, EO-gas, Radiation, or Plasma for sterilization, or any other sterilization method for this product.

**IMPORTANT:**

Prior to usage, let the product cool to room temperature. It is the sole responsibility of the user to maintain the sterile condition of the product after the sterilization process.

In case the aforementioned chemicals and machines for cleaning, disinfection, or sterilization are not available, it is the responsibility of the user to validate the procedure actually applied. Also, if a sterilization method other than that described above is applied, this deviating procedure has to be validated by the user accordingly.

**Limitation of Reprocessing**

The product lifetime is depending on wear and tear, handling, application time, damages as well as the frequency of reprocessing. Therefore, a visual inspection and function test has to be done prior to each usage. (Refer to section "Visual Inspection" and Function Test). Only an undamaged product is allowed to be reused.

**Visual Inspection and Function Test**

Prior to each usage, the complete product, especially the insulation has to be inspected for pressure points and damages. A product exhibiting damages, pressure points, or questionable conditions is not allowed to be used and has to be replaced by a new one. Prior to usage, an electrical continuity test has to be done. In case the product does not pass the electrical continuity test, it is not allowed to use the product any longer and has to be replaced by a new one. Note: If the cable conductor is broken or the electrical continuity of the cable is interrupted in any other manner it can cause arcing in the patient return system or the active circuit. This can lead to the patient suffering burns or the outbreak of a fire. If damaged, do not use these products. Never lay this product on the patient or in their direct vicinity.

**Exclusion of Repair and Modification**

Unauthorized repairs are strictly prohibited.

The product must not be bent or manipulated in any other way.

This may result in serious injury to the patient, user or third parties.

Unauthorized modifications and repair work are strictly prohibited and will result in the immediate loss of the manufacturer's warranty.

## Packaging, Storage and Transport

Store in a clean and dry environment. Storage only in protective cases with individual compartments or individually sealed in foil. Handle with the utmost care during transport, cleaning, disinfection, maintenance, sterilization, and storage (see also section "after the application"). This especially applies to fine tips and other sensitive areas. Do not store or transport the product together with sharp-edged or pointed objects. Maintenance of the sterile condition after the sterilization process is the sole responsibility of the user.

## Warranty

Our products meet the highest quality standards.

Liability for products that have been modified, that have not been used in accordance with their intended use, or that have been improperly treated or handled, or in the event of a deviation from these instructions for use, is excluded.

Furthermore, DEWIMED Medizintechnik GmbH does not assume any liability for intentional, accidental damage or damage resulting from use or handling.

In case of repairs on our products by companies, which are not authorized by us, any warranty is void.

Unauthorized repairs are strictly prohibited.

## Returns

Returned products only are accepted, if they are marked as "hygienically safe" and have been packed safely for the transport. Use our return merchandise authorization & decontamination certificate for returns.

## Disposal

Disposal of the product, their packaging material as well as any accessory has to be done according to the applicable country specific requirements, regulations and laws.

Additionally, the applicable requirements of the respective clinical place(s) in regards to disposal of medical devices have to be followed as well.

## Regulatory Remark


Report serious incidents with the product to us. If you are a user in the European Union, also report incidents to the responsible authority in your Member State.

## About these Instructions for Use

Throughout the period of usage these IFU must be kept freely accessible for the user. For a current revision of this IFU, please contact our customer service.

Changes reserved.

## Symbols used on the label

	Article number		Lot-number		Follow instruction for use
	Medical device		Manufacturer		CE-Mark
	Quantity in packaging		Manufacturing date		Non-sterile
	Prescriptive medic device (USA)				

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