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
IFU-39-V002 | 28.04.2023

## INSTRUCTIONS FOR USE STERILIZATION CONTAINER SYSTEM

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Tel. +49(0)7461-92393-0  
Fax. +49(0)7461-92393-33  
www.dewimed.de

#### Products

These instructions for use apply to the DEWIMED sterilization container system of risk class 1: Article no.: 88.XXX.XX

**Important information**  
The instructions for the correct use of Sterilization containers are described below.  
It is recommended to read the instruction Manual prior to the usage in order to avoid the risks, as much as possible, that may occur from improper usage of the container.

#### Intended Use

Sterilization container systems are systems that provide sterile storage of surgical hand tools. The system, which is suitable for reuse, is an effective packaging method that ensures sterilization, storage and transportation of surgical and orthopedic instruments. Various sizes and configurations are available. The sterilization container system is designed in accordance with Vacuum Steam Sterilization (If different sterilization method is to be applied, please contact the manufacturer or authorized distributor). The sterilization container system provides the opportunity to store the instruments, devices, etc., in a sterile form, until they are used again after the sterilization process. According to the sterility and stability test results, if the lid of sterilization container systems is not opened, they are kept sterile for 6 months

#### Sterilization container models

##### Standard-Model

This model is suitable for usage with filters;

- Paper Filter
- Textile Filter
- PTFE (Teflon) Filter
- PP (Polypropylene) Filter

##### Window Model


The Window Model Sterilization container provides the opportunity to control the condition of both product and sterilization via the transparent PPSU window system on the lid. The requirements of the sterilization process are provided by the filter, valve system or B3 in the box. Since this model has no perforation on the lid, it should be used with perforated boxes.

##### Bio-Barrier-Model

This system provides the sterilization to be fulfilled without using any filters. During the Vacuum Steam Sterilization Process, the dry air is discharged from the container via the opening of the valves through the pre-vacuum. In the sterilization stage, the valves are opened in the opposite direction and the infed of the hot air vapour is provided. After sterilization, the hot air is released again by vacuum and the valves are balanced so that there is no air flow. Thus the inner ambient of the container keeps its sterility.

##### Modell B3

The system, which consists of a combination of 3 different modules in one cover, allows to be used as PTFE Filter, Valve or Bio-Stop (Pasteur) system with the preferred module.

#### Maintenance

Proper cleaning methods shall be used for the cleaning of the sterilization container system prior to the usage. Sterilization container systems are suitable for automatic washing. Neutral Enzymatic solutions should be used during the washing process. The preferred filters for filtered models, shall be checked after each cleaning and before each sterilization. In case of any kind of damages on the filters such as cracks, scratches etc. or if the filters are wet etc. the filters shall never be used.

**Warning:** Sterilization Container Systems shall only be used with the components that are approved by our company. Our company shall not be responsible for the damages that may occur due to the usage of components unapproved or inappropriate. Ensure that all components of the Sterilization Container System work properly before each usage.

#### Unacceptable situations

- Damaged lock mechanism
  - Deformation or rupture of the internal lid gasket
- Use given values in the table for the replacement of the gasket.

#### Control and replacement of the gasket

- Clean the lid groove with the help of acetone in order not to have any gasket residues in the groove.

- Apply the silicone adhesive (Henkel-Sista® is recommended) to the gasket groove, which is cleaned properly.
- Cut the gasket related with the dimensions given in Table-1 and place it into the groove (flat surface facing up) and wait 12 hours to have the gasket dried.

Lid	X long	X Large	(1/1)	(3/4)	Fla/ (1/2)	Dental	Dental 1/2	Large Dental	Mini	Endo-Modell
Gasket Dimensions	185 cm	180 cm	158 cm	136 cm	101 cm	89 cm	72 cm	118 cm	77 cm	130 cm

- Check the gaskets in the filter holders.
- Ensure that the filter holders have been fitted exactly.
- The filter holders must fit via the "click" sound.
- Do not use the damaged, torn or wet filters.
- Check the flexibility of valve springs in B3 Model Valve System.
- Check the flexibility of the valves by pushing the valve shaft on the module by pressing with finger.
- In the Bio-Barrier Model, control the spring of the valve mechanism with a thin-tipped material.
- Do not use the mechanisms that have lost their flexibility.

The model conveniences of the filters and their usage periods are given in the table below:

Model	PP Filter <sup>2</sup>	PTFE Filter	Textil Filter	Paper Filter
Standard-Model	✓	✓	✓	✓
Window-Model	✓ <sup>1</sup>	✓ <sup>1</sup>	✓ <sup>1</sup>	✓ <sup>1</sup>
Bio-Barrier-Model	-	-	-	-
Model B3	✓	✓	-	-
Plasma-Model	✓	-	-	-
Dental-Model	✓	✓	-	✓
½-Dental-Model	✓	✓	✓	✓
Endo-Model	✓	✓	✓	✓
Mini-Model	✓	✓	✓	✓
Flat-Model	✓	✓	✓	✓
Filter by time	Single use <sup>3</sup>	2000 SC <sup>4</sup>	50 SC <sup>4</sup>	Single use <sup>3</sup>

<sup>1</sup> The filters are placed in the box in Window Model System

<sup>2</sup> PP: Polypropylene Filter

<sup>3</sup> Single-Use: Needs to be replaced after each use.

<sup>4</sup> Sterilization Cycle: It increases by 1 after each sterilization process.

SC: Sterilization Cycle

#### Cleaning and disinfection

Incorrect cleaning and disinfection may cause corrosion and erosion on the surface.

**Warning:** Use alkaline-free, caustic-free (NaOH) and carbonate-free neutral enzymatic detergents or washing. The final rinsing water to be used before sterilization shall be de-ionised (in accordance with EN 285 Appendix –B). Before washing, remove all the removable parts of the sterilization container. Remove the filter holder and the filter for the filtered model; filter module, Bio Stop Module and valve module for B3 Model. Remove the paper labels and aluminium labels.

#### Manual washing

Wash using cold water for at least 2 minutes. While rinsing use a soft brush and a lint-free cloth to ensure that the sterilization container is cleaned from dirt and dust etc. Immerse the container in neutral enzymatic detergent and soak for 5 minutes at least. Then apply the detergent on the container with a soft brush and a lint-free cloth. Bio-Barrier model is not suitable for manual washing.

**Warning:** In order to use the neutral enzymatic cleaning solutions properly, it is recommended to follow-up the manufacturer's instructions while preparing the solution (Proper Temperature, Concentration etc.).

#### Mechanical Washing

Mechanical washing of Sterilization Container is only recommended if the disinfectors has a special program for the aluminium boxes. The used program of the washers and the disinfectors for stainless steel surgical instruments, Cannot be applied to the Sterilization Containers.

**Warning:** Do not wash the container with Alkaline and Acidic solutions in automatic washers.

Process Sequence	Time (min)/ Temperature	Water Quality	Detergent Features
Pre-Wash Rinse	1/ Cold	Tap water	n/a
Wash 1 (Manual)	2/ Cold	Tap water	Neutral pH 7
Wash 2 (Mechanic)	5/ 40°C	Tap water	Neutral pH 7 Enzymatic
Final Rinse	2/ 40°C	Distilled Water	n/a
Drying	40 min/ 93°C	n/a	n/a

**Warning:** Distilled water should be used in the final rinsing process of the Sterilization Containers due to the fact that the mineral salts in mains-water may cause spots and erosion over time during sterilization.

- Remove the Sterilization Container's box from the lid prior to washing.
- Place the box vertically and upside down into the disinfectant in order to prevent the water accumulation.
- The inside of the container lid should be positioned at the bottom of the machine with the locking mechanism folded into the lid.
- Remove all the components of the container and clean them separately by placing into the basket (if applicable).

#### Loading the equipments

- Place the equipments/materials to be sterilized in the container by leaving at least 2 inches ( approx. 5 cm) gap between the lid holes and the load so that these don't prevent the function of the filters.
- If the perforations are clogged while the container is being loaded before sterilization, the vacuum in the pre-vacuum steam sterilization causes the deformation of the container.
- It is recommended to use a basket when loading the container. Baskets are used for the grouping and stacking of surgical and orthopaedic equipments (organising) in Sterilization container and also for washing these equipments.

Loading should be done according to the values given in the table below. When the values in the table are exceeded, the materials are likely to remain wet. The Sterilization process will not be complete.

Material	Container Type	Maximal Total Weight <sup>1</sup>	Relevant Standard
Surgical Instrument	1/1	10 kg	EN 868-8
Surgical Instrument	3/4	7 kg	EN 868-8
Surgical Instrument	1/2	5 kg	EN 868-8
Surgical Instrument	Mini, Dental & Special Models	3 kg	EN 868-8
Textile	All models	8kg <sup>2</sup>	EN 285

<sup>1</sup> Maximum total weight : Container weight + basket or drape weight + Material weight.

<sup>2</sup> Ensure that the container perforations do not close when the textile materials are loaded. When light pressure is applied on it, it should be positioned so as to maintain its flexural stiffness. Otherwise, the steam cannot penetrate into the container and there won't be any sterilization.

#### Sterilization

- Ensure that the lid is fully closed after the equipments are loaded. Ensure that the lock system works correctly.
- Place the lid parallel on the box and lock the lid. Otherwise, the box will not be locked properly, thus there will be an unbalanced force to be applied and the locking system won't be functioning. As the result of this there shall not be a proper sterilization.

**Warning:** Do not wrap the sterilization container in any kind of packages. In case the perforations are clogged during the sterilization, vapour and steam won't be discharged. Thus the container will be damaged and there shall not be a proper sterilization.


**Warning:** It is recommended to use either security seal or automatic lock system to ensure the sterility.

Sterilization Type	Time 121°C	Time 134°C	Drying Time
Gravity Sterilization	30 Minutes	15 Minutes	15-30 Minutes
Pre-Vacuum Steam Sterilization	-	5 Minutes	15-20 Minutes

- Prevent the products to face big temperature differences in the environment after sterilization process.
- Wait until the container cools down to the room temperature before unloading the autoclave. In case of sudden temperature changes sweating on interior surfaces of the container will occur.

- If the loading weights given in Table 4 are not followed, the sterilizer may not fully perform the drying function. In this case the equipment and materials will not be considered as sterile.
- Place the containers having larger volumes to the bottom and the ones with smaller volumes to the top, in case more than one container is to be put into the sterilizer. Do not exceed the maximum allowed height of 45cm.

**⚠ If the equipments are still humid or wet although all of these instructions are followed, contact immediately with the manufacturer of the sterilizer machine.**

#### Transportation

- Always hold from the handles of the container during loading, unloading and transportation.
- Do not carry the container holding from the lock mechanism.

**⚠ Carrying the container holding from the locking mechanism may cause a defect on the locking system and as the result of this there may be air leakage, which will lead to the loss of sterilization.**

Safe carrying weights of Sterilization Container System is given the table below.

Size	Height	Safe Carrying Weight	Stack Capacity (Up to 45 cm)
1/1 - 3/4 - 1/2	100 mm	20 kg	4 Stück
XL	120 mm	22 kg	4 Stück
1/1 - 3/4 - 1/2	135 mm	23,5 kg	4 Stück
1/1 - 3/4 - 1/2	150 mm	25 kg	4 Stück
XL	180 mm	28 kg	4 Stück
1/1 - 1/2	200 mm	30 kg	4 Stück
1/1 - 1/2	260 mm	36 kg	4 Stück

#### Storage

Recommended proper storage conditions of Sterilization container system are as follows.

- Temperature: (-20) - 60° C
- Humidity: % 10 - 70 Rh
- Pressure: Normal Atmospheric Pressure (Refer to: DIN 58953-9)

**⚠ Storage times of sterilization container systems depend on ambient conditions. Storage conditions and times must be approved by qualified personnel.**

#### Troubleshooting table

Fault	Reason	Solution
Sterilization Excessive moisture from inside the container.	The temperature of the materials before sterilization is too low	Sterilize materials at 20°C sterilization at minimum initial room Temperature
	The drapes are very damp	The drapes to be placed in the sterilization container should be dry
	Containers to be sterilized are not properly placed in the sterilization unit.	The total height should be 45 cm, with the heavier containers at the box. It should be preferred that the polymer lid container are at the top.
	Improper cooling after sterilization.	Select the appropriate sterilization parameters.
	Autoclave specifications do not meet the DIN EN 285 standard	When using the autoclave, consider the manufacturer's recommendations. Check the drying time and the steam quality. Calibrate if necessary (See Table 4)
The sterile container is overloaded		Check the appropriate weight.
	The idle loop and vacuum test are not performed daily before starting sterilization	Perform idle cycle and vacuum tests before sterilization. Make sure your sterilizer is working properly.
Locking the lid failure	Locks have been used instead of handles to carry the sterilization	Please contact Dewimed Medizintechnik GmbH. Container with a broken lock is not suitable for sterilization.

	container system	
	The sterilization container system has had a physical impact on the lid locks	
B3 Biostop stand parts are not mounted on the cover	B3 Biostop components or lid deformed.	Please contact Dewimed Medizintechnik GmbH
During the sterilization process, the B3 Bio-stop barrier components were disassembled from the chamber.	B3 Bio-stop components are not properly locked.	Please contact Dewimed Medizintechnik GmbH

#### More information

Only validated processes may be used for the reprocessing of medical devices

#### Explanation of symbols

<b>REF</b>	Article number	<b>MD</b>	Medical device	<b>QTY</b>	Quantity in packaging
<b>LOT</b>	Lot-number		Manufacturer		Manufacturing date
	Follow instruction for use	<b>CE</b>	CE-Mark		Non-sterile
	Attention	<b>R<sub>only</sub></b>	Prescriptive medical device (USA)		

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