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Monopolare HF-Elektroden

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IMPORTANT INFORMATION PLEASE READ BEFORE USE

Scope

Non-sterile, reusable product

Monopolar HF Electrodes Art. No.: 90.010.12 to 90.025.54, 90.026.02 to 90.026.88, 90.027.01 to 90.027.66, 90.030.00 to 90.030.15

Maximum rated voltage of accessory (Umax):

Article No.	Umax
90.025.32 to 90.025.54; 90.027.01 to 90.027.12	0,5 kVp
90.026.02 to 90.026.04	1,3 kVp

Article No.	Umax
90.010.12 to 90.025.30; 90.026.05 to 90.026.88; 90.027.13 to 90.027.66; 90.030.00 bis 90.030.15	4,3 kVp

See also labeling or catalogue.

In any combination with another electrosurgical accessory, the maximum rated voltage of the combination corresponds to the lowest rated voltage of the accessories used. (See also section "Intended Use") See these instructions for Use, the label, or the current product catalog for the maximum rated voltage of the product. In case of uncertainties, contact the manufacturer. Prior to usage, read the entire IFU of this product and of any used accessory as well as HF-generator and HF-neutral electrode (monopolar application). All requirements, safety notices and warnings included in the respective IFU have to be followed strictly. This medium is not intended for users in the USA.

Intended Use

The product is intended for open or endoscopic surgery and is used for cutting and coagulating biological tissue.

Only for use by skilled medical professionals.

Indication:

Product intended for use in open or rather endoscopic surgery.

Contraindication:

Product is not intended for direct contact with the heart or the central circulatory or central nervous system. Application of high-frequency current may interfere with cardiac pacemakers and in vivo heart defibrillators, so affected patients must consult a cardiologist prior to the intervention.

Connection and Activation:

Prior to application of monopolar products, it has to be safeguarded that an HF neutral electrode has been applied correctly to the patient and that this device is connected properly to the respective HF-generator. The product is connected by a suitable HF-cable or a suitable HF-handle to the respective exit of an HF-generator. Activation is done using the buttons of the HF-handle or with the foot switch of the HF-generator.

Combination/Compatibility:

Prior to use, compatibility of the product with the intended HF-handle/HF-cable and HF-generator has to be verified. The frequency of the used HF-generator shall not exceed the maximum frequency of 4 MHz and not exceed the maximum rated voltage of the accessory (see section "Scope"). In case of uncertainties, contact the manufacturer of this product or the manufacturer of the used accessory or HF-generator.

The monopolar HF electrodes can be operated on generators from the manufacturers listed below:

ERBE	KLS Martin	EMED	BOWA

Follow the instructions, safety notices and warnings included in the User Manuals or IFU of the used accessories and used HF-generator.

It is recommended to use a smoke evacuation system.

IMPORTANT: Handle with utmost care. This does not only apply for the duration of the surgery but also for the complete duration of storage, processing and transport as well as during the process of connecting the product with the HF-accessory and HF-generator. This applies especially for the thin components and other sensitive areas, e.g. the insulation. **Prior to first use** as well as **after each use**, the product has to be cleaned, disinfected and sterilized (see section "(Re-) Processing: Cleaning, Disinfection and Sterilization"). Improper use immediately will result in loss of warranty. Liability for any damages incurred will not be accepted.

Safety Notices – WARNING

See these instructions for use, the label, or the current product catalog for the maximum rated voltage of the product. Follow the aforementioned allowed frequency of HF current. Prior to each application, the product has to be cleaned, disinfected and sterilized according to a validation procedure (EN ISO 17665). Prior to usage a visual inspection and function test has to be done (see section "Visual Inspection and Function Test"). In case damages, deformation or similar is detected on the product, it is not allowed to use the device. It has to be replaced by a new one. At least one (1) cleaned, disinfected and sterilized backup product has to be available. It is the responsibility of the user to determine the appropriate product size and product type according to their professional judgement and based on the patient's specific indication, preferred surgical technique and history, etc. Prior to use, ensure that the product is firmly inserted in the HF handle/HF cable. This must be done carefully, in order to avoid damages on the product and/or injuries of patient, surgical personnel or third party. Excessive force can damage the product. Therefore, the product has to be observed during the complete application.

Exclusion:

Do not activate the product as long as it is in contact with metal objects and/or optics. During an electrosurgical intervention the patient must not come into contact with grounded metal objects such as surgical desk frames, instrument trays etc. Pay attention that no flammable substances are present in immediate vicinity as otherwise there is danger of explosion.

(Refer to section "General Safety Notices for HF-Technology").

After switching off the electrosurgical current, the product tip may still be hot and can lead to unintentional burns. Do not activate the product uninterruptedly during a longer period. It is the responsibility of the user to apply low performance settings on the HF-Generator in order to achieve the envisaged effect for the respective intervention. During the surgery, mechanical forces may lead to deformation or signs of wear of the product. If, abrasion of the product is detected during the application, the product has to be exchanged by a new one. Abrasion, adhesion of tissue, discoloration, sooting etc. do not represent a reason for complaint and do not permit claiming the manufacturer's warranty.

General Safety Notices for HF-Technology (Excerpt)

In addition to the acknowledged benefits of HF-Technology, the application includes several risks that have to be attended: Improper use and non-observance of the IFU can lead to unintentional burns of the patient, user or third party.

Continued further education of the surgical personnel is recommended.

a) Environment

Pay attention that no flammable substances (anesthetics, oxidizing gases, endogenous gases etc.) are present in immediate vicinity as otherwise there is danger of explosion. Only use non-

flammable cleaning or disinfection agents. All oxygen connections must be tight and leak-proof during the procedure.

b) Patient Positioning and Preparation

Ensure proper patient positioning, i.e., only use insulating surgical drapes that are dry, absorbent and liquid-tight. Isolate conductive surfaces and contact points towards the patient. Dry tissues of pulp are necessary to use for skin and breast folds as well as between the extremities. Prior to application, remove any liquids that potentially accumulated in body cavities. Only use non-flammable disinfectants. Do not use alcohol-based tinctures. Only use non-conductive irrigation fluid, if medically possible. Attend the requirements on irrigation fluid for monopolar and bipolar products. Prior to application, remove any body jewelry from the patient. Putting a band aid over the body jewelry is not sufficient!

c) Connections

Prior to application, make sure the product is connected properly to the HF-generator and make sure the correct power and performance setting is adjusted and displayed. It is the responsibility of the user to apply low performance settings on the HF-Generator in order to achieve the envisaged effect for the respective intervention.

d) HF-Neutral Electrode for Monopolar Application

In case of monopolar application, select a HF-neutral electrode suitable for the patient, apply it correctly and connect it properly with the respective HF-generator. Follow all instructions for proper application of the HF-neutral electrode, incl. patient protection and patient monitoring, monitoring of the HF-neutral electrode and all further provisions, safety notices and warnings included in the IFU of the HF-neutral electrode.

e) Patient Reactions

All electrosurgical devices potentially can cause muscle stimulation during the application. The design of this product minimizes the risk of this undesirable effect. Nevertheless, muscle stimulation can lead to an unexpected movement of the patient in the surgical field.

f) Handling HF-Accessory

Make sure the accessory is compatible. Do not touch the instrument tip during the complete application. As long as the product is not applied, place it on a dry, clean and non-conductive and well-visible surface, that is not in contact with the patient. Never store product on the patient. Unintended activation of the product can lead to unintentional burns or other injuries of patient, user or third party. Never wrap cables around the patient and never lay cables over the patient. Apply only short activation times. Keep on longer breaks between activation phases. Only adjust low power settings.

g) Completeness of the System

At the end of the surgery, confirm completeness of the system.

(Re) Processing: Cleaning, Disinfection and Sterilization

Maximum number of reprocessing cycles

Due to the design, materials used, intended use as well as wear and tear, a maximum limit of performable reprocessing cycles cannot be determined.

When applied according to the Intended Use, the product underlies natural wear and tear, considering manner and duration of the application as well as manner and frequency of reprocessing. Therefore, a visual inspection and function test has to be done prior to each usage. (Refer to section "Visual Inspection" and "Function Test"). Visual inspection and function test, especially the condition of the complete cable, insulation, port, plug component and pins are decisive for whether the product is allowed to be applied again.

Time requirements for cleaning and disinfection

Preparation for cleaning, pre-cleaning and automated cleaning and disinfection has to be done immediately after the application, however not later than 1 hour after the application. Avoid idle time.

After the Application

Clean and disinfect product immediately after the application. However not later than 1 h after application. After application of the product, deposit the product carefully (protection of lifetime of the product). After application, separate contaminated product and deposit it in a suitable container (deposit means "do not drop"). Immediately remove gross stain. Immediately mark damaged or defective products. Accessory, that do not fit for the sieves of the cleaning and disinfection device (CDD), shall be deposit separately in suitable containers. Close firmly all disposal-/transport containers, in order to avoid drying of stain. Organize transport of contaminated products the way that contamination of the transport ways and environment also is avoided (closed transport). Unused reusable products have to be reprocessed as well. (Refer to AKI (German organization for instrument reprocessing), Red Brochure, page 30-32). Take care, all and any transport containers are cleaned and disinfected after transport as well.

Validation of (Re)Processing

The following validated processing procedure is recommended.

Equivalent deviant processes are possible. Then, it is the sole responsibility of the user to safeguard the suitability of the actually applied procedure by suitable means (e.g. validation, routine monitoring, verification of material compatibility etc.). Automated cleaning and disinfection always are preferable.

The following procedure was validated according to EN ISO 17665. Additional applicable processing requirements specific to the respective clinical place (operator) as well as national or country specific regulations have to be followed as well. Never use sharp objects for cleaning. Disinfectants always have to be rinsed and removed carefully.

Preparation for Cleaning

Remove product from their packaging. Place it in a container provided for cleaning. It is not necessary to disassemble the product.

Pre-Cleaning

Immediately after the application is completed, pre-clean the product. This however not later than 1 hour after the end of the surgery. Use tap water (potable water quality) (<40°C) and aldehyde-free, non-fixing disinfectants if applicable.

Thoroughly remove surface stain with a soft brush or synthetic fleece, as otherwise particles or dried secretions may adhere. This could make subsequent cleaning and sterilization difficult or impossible. Ensure that areas difficult to reach are cleaned thoroughly and rinsed several times. Cavities and lumen have to be rinsed intensively using at least 3x20 ml cold tap water (<40°C) with the aid of a rinsing adapter (e.g. from the company Medisafe), or with a syringe or with a water jet pistol (> 30 Sec.). This pre-cleaning step always has to be done prior to the manual cleaning or cleaning with the cleaning and disinfection device (CDD).

Manual Cleaning and Disinfection

Prepare an immersion bath with a suitable fluid cleaning agent.

Use a cleaning agent that is compatible to the disinfectant and suitable for immersion baths. Follow the instructions of the manufacturer of cleaning agent and disinfectant. Only use agents suitable for medical devices made from metal and plastics with a pH-value between 5.5 and 12.3. Recommendation: Cleaning agent gigazyme® (Schülke & Mayr) and disinfectant Korsolex Plus.

Do not use high alkaline cleaning agents.

These will impair the lifetime of the product.

- Prepare immersion bath with cleaning agent according to the specific cleaning agent IFU.
- Prepare a separate immersion bath with disinfectant according to the specific disinfectant IFU.
- Immerse the product completely in the ultrasonic bath with cleaning agent (e.g. 0.5% gigazyme®)
- Clean product in ultrasonic bath using a sonication time of 5 Min. and a frequency of 35 kHz.
- Follow all instructions set forth in the IFU of the cleaning agent, disinfectant and ultrasonic bath.
- Ensure the product will not touch other products or parts in the ultrasonic bath.
- Ensure skin shadows in the ultrasonic bath are avoided.
- Then, clean product with a soft brush under cold running town water (<40°C).
- Intensively rinse cavities and lumen with a water jet pistol (>30 Sec.) or similar for at least 1 Min.
- Afterwards, rinse product thoroughly for at least 1 min. under running tap water (>40°C) in order to remove any residues of the cleaning agent.
- Inspect product visually on remaining stain.
- In case stain is still present, repeat aforementioned cleaning steps as long as it needs until no soil is present.
- Afterwards: immerse product completely in a disinfectant bath including e.g. Korsolex Plus 3%, for at least 15 Min.
- Follow the manufacturer's data for residence time.
- Ensure the disinfectant will contact all areas on the product.
- Rinse cavities and lumen several times, that means at least 3x with 20 ml each of disinfectant bath fluid.
- Afterwards: Rinse product thoroughly for at least 1 min. with demineralized cold water, in order to remove all disinfectant residues.
- Additionally: Rinse all narrow and areas difficult to access, all cavities and lumen with a syringe several times (at least 3x) using each time 20 ml cold demineralized water.
- Dry product with a lint free wipe and sterile compressed air.
- Dry cavities, lumen and channels with sterile compressed air.

Automated Cleaning and Disinfection

Only use cleaning and disinfection devices (CDD) with proven efficiency according to EN ISO 15883. Follow data of the manufacturer of the cleaning and disinfectant machine. Only use agents suitable for medical devices made from metal and plastics with a pH-value between 5.5 and 12.3.

Recommendation: neodisher® mediclean forte (Dr. Weigert GmbH & Co. KG).

Apply program for thermal disinfection. Follow instructions and data regarding program course and machine. Product has to be stored safely and protected against mechanical damages during automated cleaning and disinfection.

Do not clean together with sharp edged or pointed objects. Deposit product in a suitable rinsing basket. Follow data for loading of the cleaning and disinfection device (CDD). Use rinsing adapters for products with lumen and connect them according to the instructions in the User Manual of the cleaning and disinfection device (CDD).

Cleaning Program

Start the program course with following parameters:

- 1 min. pre-rinsing with cold water
- Emptying
- 3 min. pre-rinsing with cold water
- Emptying
- 5 min. cleaning at 55°C with 0.5% alkaline cleaning agent
- Emptying
- 3 min. neutralization with warm tap water (>40°C) and neutralizer (0.1% Neodisher® Z)
- Emptying
- 2 min. interim rinsing with warm demineralized water (>40°C)
- Emptying

Disinfection Program

Automated thermal disinfection considering national requirements regarding AO value (see EN ISO 15883, AO value >3000).

- 5 min. cleaning at 92°C +/-2°C

Drying

- 30 min. at 90°C
- Remove rinsing adapter
- At the end of program course, remove product and inspect it on remaining stain. In case of residues, repeat automated cleaning and disinfection step as long as it takes until stain is no longer present.

Dry cavities and not sufficiently dried areas with sterile compressed air <2 bar. Immediately after removal of product and immediately after additional drying on a clean plate, put product in a single use sterilization packing (double packing) from paper or foil or put product in a sterilization container. Respect requirements for sterilization packaging according to EN ISO 11607 and EN 868.

Sterilization

Only products that have been cleaned and disinfected are allowed to be sterilized. Only apply steam sterilization in autoclave (fractioned pre-vacuum with sufficient product drying) for this product.

Adjust sterilization parameters:

- Minimum 134°C and maximum 137°C in saturated steam.
- Holding time at least 5 Min. until max. 20 Min.
- Drying in vacuum for at least 10 Min.

- Sterilizer (Class B) according to valid national standards and regulation (e.g. EN 13060 or EN 285).

Example: Sterilizer Class B, manufacturer: Tuttnauer.

Respect data of the Sterilizer manufacturer regarding load, handling and drying times.

Exclusion:

Do not apply hot air, EO-gas, Radiation or Plasma for sterilization, or any other sterilization method for this product.

IMPORTANT:

Prior to usage, let product cool to room temperature. It is the sole responsibility of the user to maintain the sterile condition of the product after the sterilization process. In case the aforementioned chemicals and machines for cleaning, disinfection or sterilization are not available, it is the responsibility of the user to validate the procedure actually applied. Also, if a sterilization method other than described above is applied, this deviating procedure has to be validated by the user accordingly.

Limitation of Reprocessing

The product life time is depending on wear and tear, handling, application time, damages as well as frequency of reprocessing. Therefore, a visual inspection and function test has to be done prior to each usage. (Refer to section "Visual Inspection" and "Function Test"). Only an undamaged product is allowed to be reused.

Visual Inspection and Function Test

Prior each use, check the entire product, especially the insulation and product tip, for pressure marks and damage. A product exhibiting damages, pressure points or questionable condition is not allowed to be used and has to be replaced by a new one. During and after application, tissue may adhere to the product, or sooting may be present on the distal end of the active electrode.

Such adhesions or sooting do not represent a reason for complaint and the product has to be exchanged by a new one. Due to longer application time, mechanical forces or plasma seam or similar, the product may exhibit deformation or abrasion of the insulation material. Also, such aspects do not represent a reason for complaint and the product has to be replaced by a new one. Blockage of the suction channel (if applicable), does not represent a reason for complaint. A product, that exhibits a blocked suction channel (if applicable), has to be replaced by a new one. Prior to usage an electrical continuity test has to be done. In case the product does not pass the electrical continuity test, it is not allowed to use the product any longer and has to be replaced by a new one.

Exclusion of Repair and Modification

Unauthorized modification and repairs (e.g. bending) are strictly prohibited.

Especialmente, teniendo un gancho como activo part, are never allowed to be bent.

This could lead to severe injuries of patient, user or third party. Manufacturer's warranty immediately is terminated in case of modification or repair of the product.

Packaging, Storage and Transport

Store in a clean and dry environment.

Storage only in protective cases with individual compartments or individually sealed in foil. Handle with utmost care during transport, cleaning, disinfection, maintenance, sterilization and storage (see also section "after the application"). This especially applies for fine tips and other sensitive areas. Do not store or transport the product together with sharp edged or pointed objects. Maintenance of the sterile condition after sterilization process is the sole responsibility of the user.

Warranty

DEWIMED Medizintechnik GmbH only provides inspected and faultless products. Our products satisfy the highest quality standards. Liability and warranty are excluded for all products that have been modified in any way, not applied according to their Intended Use, or that have been handled or applied improperly, or in case of any other deviation from the instructions set forth in this IFU. Furthermore, DEWIMED Medizintechnik GmbH denies any liability for any accidental, intentional damage or for a damage or loss arising out of handling or application of the product. Additionally, all liability and warranty are terminated in case our product was repaired by a company that has not been authorized by us. Unauthorized repairs are strictly prohibited.

Return

Returned products only are accepted, if they are marked as "hygienically safe" and have been packed safely for the transport. Use our return merchandise authorization & decontamination certificate for returns.

Disposal

Disposal of the product, their packaging material as well as any accessory has to be done according to the applicable country specific requirements, regulations and laws. Additionally, the applicable requirements of the respective clinical place(s) in regards to disposal of medical devices have to be followed as well.

Regulatory Remark


Report serious incidents with the product to us. If you are a user in the European Union, also report incidents to the responsible authority in your Member State.


About these Instructions for Use

Throughout the period of usage these IFU must be kept freely accessible for the user. For a current revision of this IFU, please contact our customer service.

Changes reserved. This medium is not intended for users in the USA.

Explanation of symbols

	Attention		Manufacturer		Manufacturing date
	Medical device		Non-sterile		Quantity in packaging
	Article number		Lot-number		Follow instruction for use

	Keep dry		Keep away from sunlight
	CE-Mark of the Notified Body mdc medical device certification GmbH, Stuttgart, Germany		Prescriptive medical device (USA)
		IFU-20-V003 16.02.2024	
ESPAÑOL		Electrodos HF monopolares	
INSTRUCCIONES DE USO		Electrodos HF monopolares	
	Dewimed Medizintechnik GmbH Unter Hasslen 14 78532 Tuttingen ALEMANIA	Tel. +49(0)7461-92393-0 Fax. +49(0)7461-92393-33 www.dewimed.de	
			
			
POR FAVOR, LEA LA INFORMACIÓN IMPORTANTE ANTES DE USAR			

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Marque CE de l'organisme notifié mdc medical device certification GmbH, Stuttgart, Allemagne