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HF-Kabel/ Adapter Art.-Nr.: 90.100.40, 90.100.48, 90.190.30, 90.190.50, 90.192.30, 90.192.50, 90.193.30, 90.193.50, 90.196.30, 90.196.50, 90.197.30, 90.197.50, 90.199.30, 90.199.50, 90.201.50, 90.202.50, 90.204.50, 90.205.30, 90.205.50, 90.206.30, 90.206.50, 90.208.30, 90.208.50, 90.209.50, 90.210.50, 90.212.50, 90.213.30, 90.213.50, 90.214.30, 90.214.50, 90.215.30, 90.215.50, 90.220.30, 90.220.50, 90.221.30, 90.221.50, 90.222.30, 90.222.50, 90.223.30, 90.223.50, 90.224.30, 90.224.50, 90.225.30, 90.225.50, 90.226.40, 90.226.50, 90.227.40, 90.227.50, 90.229.30, 90.229.50.

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# Instruction for use Electrosurgical Cables and Adapters

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Copyright by DEWIMED Medizintechnik GmbH  
All rights reserved  
Subject to technical changes without notice

Revision: 03

Date: 2016-05-30

## ENGLISH

### 1 Scope

Electrosurgical cables/ adapters Art. No.: .90.100.40, 90.100.48, 90.190.30, 90.190.50, 90.192.30, 90.192.50, 90.193.30, 90.193.50, 90.196.30, 90.196.50, 90.197.30, 90.197.50, 90.199.30, 90.199.50, 90.201.50, 90.202.50, 90.204.50, 90.205.30, 90.205.50, 90.206.30, 90.206.50, 90.208.30, 90.208.50, 90.209.50, 90.210.50, 90.212.50, 90.213.30, 90.213.50, 90.214.30, 90.214.50, 90.215.30, 90.215.50, 90.220.30, 90.220.50, 90.221.30, 90.221.50, 90.222.30, 90.222.50, 90.223.30, 90.223.50, 90.224.30, 90.224.50, 90.225.30, 90.225.50, 90.226.40, 90.226.50, 90.227.40, 90.227.50, 90.229.30, 90.229.50.

Maximum rated voltage of accessory:

Designation	U <sub>max</sub>
Monopolar cable/ adapter	4,3 kVp
Bipolar cable/ adapter	0,8 kVp
Neutral-electrode-cable	4,3 kVp

See also label or catalog.

In any combination with another electrosurgical accessory, the maximum rated voltage of the combination corresponds to the lowest rated voltage of the accessories used.

### 2 Authorized use

Only skilled medical personnel are permitted to use the electrosurgical cables/ adapters.

Electrosurgical cables/ adapters are used to connect electrosurgical-accessories for monopolar or bipolar application to an electrosurgical generator.

Follow the instructions for use of the corresponding electrosurgical accessory and of the generators used.

### 3 Safety notice

#### WARNING !

See these Instructions for Use, the label or the current product catalog for the maximum rated voltage of the accessory.

If anything is unclear, contact the manufacturer. Before initial use, new electrosurgical cables/ adapters must be cleaned and also sterilized according to a validated procedure (DIN EN 17665) (In this regard see Section 4 "Cleaning and sterilization"). A function test and a visual inspection must be performed before each use (see Section 5 "Visual inspection and function test"). It must be ensured that the electrosurgical cable/ adapter is connected firmly and correctly to the generator. Throughout the complete procedure, care must be taken that no flammable substances are present in the immediate vicinity, since otherwise a danger of explosion exists.

The high-frequency-current used in electrosurgery may interfere with cardiac pacemakers and implanted heart defibrillators, and so affected patients must consult a cardiologist prior to the operation.

### 4 Cleaning and sterilization

In view of the design, the materials used and the intended use, a maximum limit cannot be defined for

the number of cleaning and sterilization cycles that may be possible.

#### Preparation for cleaning:

Remove the electrosurgical cables/ adapters from their packaging. Place them in a container/ device provided for cleaning/ sterilization.

It is not necessary to dismantle the electrosurgical cables/ adapters.

#### Manual Pre-cleaning:

The electrosurgical cables/ adapters must be disinfected immediately after each use. The electrosurgical cables/ adapters should be thoroughly cleaned with a soft brush or synthetic fleece pad then rinsed, since otherwise particles or dried secretions may adhere to them. This could make subsequent cleaning and sterilization difficult or no longer possible. It must be ensured that areas with difficult access are thoroughly cleaned then rinsed several times. Highly alkaline cleaning agents with pH levels above 10 and below 13 have no detrimental influence on the useful life. This manual pre-cleaning has to be done before the automated cleaning.

#### Automated cleaning:

Only washing machines/ disinfectors with efficiency tested according to DIN EN ISO 15883 may be used. For cleaning, use only commercially available cleaners, which are suitable for cleaning of plastics parts. Thermal disinfection must be performed for at least 5 minutes at a minimum of 90°C [194°F] to a maximum to 95°C [203°F]. The manufacturer's instructions about the dosage of cleaning agent and selection of program must be followed. If chemical disinfection has been applied, a final rinse with distilled or demineralized water must be performed. The electrosurgical cables/ adapters must be mounted securely and protected against mechanical damage during machine cleaning/ disinfection.

#### Sterilization:

The electrosurgical cables/ adapters are designed exclusively for steam sterilization in autoclaves. The electrosurgical-cables/ adapters must be sterilized at a minimum of 134°C [273°F] and a maximum of 138°C [280°F] in saturated steam during a holding time of at least 5 minutes to at most 20 minutes, then dried in vacuum for at least 10 minutes. Sterilization must be performed in accordance with DIN EN ISO 17665 (Sterilization of Medical Devices in Moist Heat). It is not permissible to sterilize the electrosurgical cables/ adapters with hot air, EO gas, gamma radiation or plasma.

**NOTE:** Before use, the electrosurgical cables/ adapters must be cooled to room temperature.

#### Limitation of reconditioning:

The useful life of the product will depend on wear, damage and frequency of reconditioning.

### 5 Visual inspection and functional test

Before each use, the insulation of the electrosurgical cables/ adapters must be inspected for pressure points or damage. In addition, an electrical continuity test must be performed before each use. It is not permissible to use damaged electrosurgical cables/ adapters or electrosurgical cables/ adapters exhibiting pressure points.

### 6 Repair and modification

It is not permissible to repair defective electrosurgical-cables/ adapters. They must be replaced by new electrosurgical-cables/ adapters. Unauthorized modifications and repairs are strictly prohibited and will entail invalidation of the manufacturer's warranty.

### 7 Packaging, storage, transportation, handling

The electrosurgical-cables/ adapters must be stored in a clean and dry environment. They should be individually stored in a protective container with individual compartments or heat-sealed in film. The electrosurgical-cables/ adapters must always be handled with the utmost care during transportation, cleaning, upkeep, sterilization and storage. This applies in particular for fine tips and other sensitive areas. It is the operator's responsibility to ensure that the sterile condition is preserved after the sterilization process.

### 8 Returns

Returns will be accepted only if they are marked as "hygienically safe" or "not contaminated" and have been securely packaged for shipping.

### 9 Disposal

The electrosurgical-cables/ adapters, the packaging material and the accessories must be disposed of in accordance with the regulations and laws specific to the country in which they are used.

### 10 About these Instruction for Use

Throughout the period of use of the electrosurgical-cables/ adapters, the Instructions for Use must be kept freely accessible for every user.