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Instructions for use Electrosurgical Forceps

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Copyright by DEWIMED Medizintechnik GmbH
All rights reserved
Subject to technical changes without notice

IFU-06-BB-GER-ENG

last update:
2017-05-23

ENGLISH

1 Scope

Bipolar electrosurgical forceps Art. No.: 90.116.11 to 90.124.18N, 90.126.16 to 90.127.22N, 90.144.16 to 90.156.19N, 90.159.18 to 90.163.19N, 90.165.19 to 90.166.22N.

Maximum rated voltage of accessory:

Designation	Article no.	U _{max}
Bipolar forceps	90.116.11 to 90.124.18N, 90.126.16 to 90.127.22N, 90.144.16 to 90.156.19N, 90.159.18 to 90.163.19N, 90.165.19 to 90.166.22N.	0,5 kVp

See also label or catalog.

In any combination with another electrosurgical accessory, the maximum rated voltage of the combination corresponds to the lowest rated voltage of the accessories used.

2 Intended use

Only skilled medical personnel are permitted to use the electrosurgical forceps.

The electrosurgical forceps are intended for open surgery and are used for grasping, manipulating and coagulating biological tissue. They are not intended for direct contact with the heart or central circulatory system. The forceps are connected by a bipolar electrosurgical cable with European flat connector to the bipolar output of the electrosurgical generator. The maximum frequency of the generator used may not exceed 4 MHz and the maximum rated voltage (see Section 1) has to be adhered to. The electrosurgical forceps are activated by means of a foot-operated switch. It is recommended that a smoke evacuation system be used.

3 Safety notice

WARNING !

See these Instructions for Use, the label or the current product catalog for the maximum rated voltage of the electrosurgical forceps.

If anything is unclear, contact the manufacturer.

Before initial use, new electrosurgical forceps must be cleaned, disinfected and also sterilized according to a validated procedure (DIN EN ISO 17665) (In this regard see Section 4 "Cleaning, disinfection and sterilization"). A visual inspection must be performed before each use (see Section 5 "Visual inspection"). It must be ensured that the electrosurgical forceps are inserted firmly in the electrosurgical-cable. This must be done carefully, in order to avoid damage to the electrosurgical forceps and/or injuries to the patient or surgical personnel. Tips may be damaged if excessive force is applied. It is not permissible to activate the electrosurgical-forceps as long as they are in contact with metal objects and/ or optics. Throughout the complete procedure, care must be taken that no flammable substances are present in the immediate vicinity, since otherwise a danger of explosion exists. After the electrosurgical current has been turned off, the instrument tip may still be hot enough to cause burns.

The high-frequency current used in electrosurgery may interfere with cardiac pacemakers and implanted heart defibrillators, and so affected patients must consult a cardiologist prior to the operation.

4 Cleaning, disinfection and sterilization

DEWIMED Medizintechnik GmbH recommends the validated treatment processes described in the following. Equivalent deviant processes are possible. It is the sole responsibility of the user to safeguard the suitability of the actual applied procedure by suitable means, f.e. validation, routine examination, verification of material compatibility, etc.

In view of the design, the materials used and the intended use, a maximum limit cannot be defined for the number of cleaning and sterilization cycles that may be possible.

Preparation for cleaning:

Remove the electrosurgical forceps from their packaging. Place them in a container/ device provided for cleaning/ sterilization.

It is not necessary to dismantle the electrosurgical forceps.

Manual pre-cleaning:

The electrosurgical forceps must be disinfected immediately after each use. For pre cleaning use water and aldehyde-free, non-fixing disinfectant if necessary. The electrosurgical forceps should be thoroughly cleaned with a soft brush or synthetic fleece pad then rinsed, since otherwise particles or dried secretions may adhere to them. This could make subsequent cleaning and sterilization difficult or impossible. It must be ensured that areas with difficult access are thoroughly cleaned then rinsed several times.

Cavities and und lumen (f. e. rinsing/suction channel) have to be rinsed intensively using at least 3 x 20 ml cold mains water (<40°C) with the aids of a rinsing adapter (f.e. from the company Medisafe), or with a syringe or with a water pressure pistol (>30 Sec.).

This pre-cleaning step has to be done prior to the further manual cleaning or prior to cleaning with the cleaning and disinfection machine.

A) Manual cleaning and disinfection:

For cleaning please prepare an immersion bath with a suitable fluid cleaning agent. Use a disinfection agent that is compatible to the cleaning agent. Take care to use a disinfection agent that also can be used for immersion baths. Follow the instructions and recommendations of the manufacturer of the cleaning and disinfection agent. Only use agents that are suitable to clean medical devices made from metal and plastics with a pH-value of between 5.5 and 12.3. It is recommended to use the cleaning agent gigazyme® (Schülke & Mayr).

1. According to the instructions of the manufacturer please prepare two separate immersion baths, one for cleaning and one for disinfection using the appropriate and respective agents.

2. Put the electrosurgical forceps completely into the cleaner (f.e. 0.5% gigazyme®) in the ultrasonic bath applying an ultrasonic time of 5 Min. and a frequency of 35 kHz. Please follow the instructions of the manufacturer of the cleaning agent. Take care that the instruments do not touch any other parts inside the ultrasonic bath. Take care that sonic shadows are avoided inside the ultrasonic bath.

3. Then, clean the electrosurgical forceps with a soft brush under cold fluent town water (<40°C). Cavities and lumen have to be rinsed intensively by using a water pressure pistol (>30 Sec.) or similar. Afterwards rinse the electrosurgical forceps thoroughly for at least 1 minute under mains water (<40°C) in order to remove any residues of the cleaning agent.

4. Inspect the electrosurgical forceps visually regarding remaining stain. In case soiling/stain is still detectable visually, repeat the aforementioned steps as long as it needs to remove any and all stain.

5. Afterwards put the electrosurgical forceps into the disinfection bath completely. Keep on the contact time set forth by the manufacturer of the disinfection agent, f.e. Korsolox Plus, 3% for 15 Min. Make sure that all areas and parts of the electrosurgical forceps are covered by the disinfectant. Cavities and lumen have to be rinsed with the disinfectant for several times using a syringe f.e. and at least min. 3 times using 20 ml.

6. Then, rinse the instruments thoroughly for at least 1 minute with demineralized water in order to remove any and all residues of the disinfectant. Regarding electrosurgical forceps that have a rinsing/suction channel, please rinse this channel at least 3 times using each 20 ml demineralized water. Rinse the channel with a rinsing adapter (f.e. from the company Medisafe) or with a syringe.

7. Dry the electrosurgical forceps with a lint free wipe and sterile compressed air. Dry the cavities and channels with sterile compressed air as well.

B) Automated cleaning and disinfection:

Only apply cleaning and disinfecting machines that are efficiency tested according to DIN EN ISO 15883. Please refer to the information of the manufacturer of the cleaning or disinfecting agent, and only use agents that are suitable to be applied for medical devices made from metal and steel and have a pH-value between 5.5 and 12.3. It is recommended to use neodisher® mediclean forte (Dr. Weigert GmbH & Co. KG). Select the program for thermal disinfection. Regarding the program course, follow the instructions of use of the manufacturer of the agents. Do not clean electrosurgical forceps together with sharp-edged or pointed objects.

Put the electrosurgical forceps in a suitable rinsing device. Connect the instruments to the rinsing device using a rinsing adapter (f.e. from the company Medisafe).

- Pre rinse with cold water for 1 Min.
- Emptying
- Pre rinse with cold water for 3 Min.
- Emptying
- Wash at 55°C for 5 Min. with a 0.5% alkaline cleaner
- Emptying
- Neutralisation with warm mains water(>40°C) and Neutralisator (0.1% Neodisher®Z) for 3 Min.
- Emptying
- Intermediate flusing: rinse with warm demineralized water (>40°C) for 2 Min.
- Emptying

Disinfection:

Proceed with the automated Thermic Disinfection according to the national requirements regarding the AO-value (see ISO 15883, AO-value >3000): 5 Min. for 92°C +/- 2°C.

Drying:

- 30 Min. at 90°C
- Remove rinsing adapter

At the end of the program course, remove the electrosurgical forceps and verify if there are any staining residues. If staining/bonded residues are

present, repeat the automated cleaning and disinfecting process as long as no visible residues are present. Dry any cavities and insufficiently dried areas with sterile compressed air < 2 bar. After the electrosurgical forceps are taken out (after additional re-drying on a clean place), immediately pack them in a single-use sterilization package (double package) made of paper/foil or put them in a sterilization container (according to DIN EN ISO 11607 and DIN EN 868).

Sterilization:

It is only allowed to sterilize instruments that have been cleaned and disinfected beforehand. Electrosurgical forceps are designed exclusively for steam sterilization in autoclaves (fractionated fore-vacuum procedure with sufficient product drying). Electrosurgical forceps have to be sterilized at a minimum of 134°C and maximum of 137°C in saturated steam during a holding time of at least 3 minutes to at most 20 minutes, then dried in a vacuum and sterilized for at least 10 minutes. Sterilizer according to valid national Norms and Guidelines (f.e. DIN EN 13060 or DIN EN 285). The sterilization process has been validated according to DIN EN ISO 17665. Please follow the recommendations and instructions of the sterilizer manufacturer regarding loading, handling and drying times. It is not permissible to sterilize the electrosurgical forceps with hot air, EO gas, gamma radiation or plasma.

NOTE: Before use, electrosurgical forceps must be cooled to room temperature. It is the sole responsibility of the user to safeguard and maintain the sterile condition of the electrosurgical forceps after the sterilization process. In case the aforementioned described and recommended chemicals and machines for manual pre-cleaning and/or automated cleaning and disinfection are not available, it is the responsibility of the user to validate their procedure. Also, if another sterilization procedure is chosen, the procedure deviant from the procedure described in paragraph 4 has to be validated by the user accordingly.

Limitation of reconditioning:

In view of the design, the materials used and the intended use of the product, a maximum limit cannot be defined for the number of cleaning, disinfection or sterilization cycles that may be possible. The lifetime of the product will depend on natural wear and tear, damage and frequency of reconditioning.

5 Visual inspection

Before each use, the insulation of the electrosurgical forceps must be inspected for pressure points or damage. It is not permissible to use electrosurgical forceps exhibiting damage or pressure points. The tip of the bipolar forceps can be tarnishing. This is not a functional impairment or a risk for the patient. The original luster can be restored by rubbing with a silver cleaning cloth prior to the cleaning, disinfection and sterilization process.

6 Repair and modification

It is not permissible to repair defective electrosurgical forceps. They must be replaced by new electrosurgical forceps. Unauthorized modifications and repairs are strictly prohibited and will entail invalidation of the manufacturer's warranty.

7 Packaging, storage, transportation, handling

The electrosurgical forceps must be stored in a clean and dry environment. They should be individually stored in a protective container with individual compartments or heat-sealed in film. The electrosurgical forceps must always be handled with the utmost care during transportation, cleaning, upkeep, sterilization and storage. This applies in particular for fine tips and other sensitive areas. It is the operator's responsibility to ensure that the sterile condition is preserved after the sterilization process.

8 Returns

Returns will be accepted only if they are marked as "hygienically safe" or "not contaminated" and have been securely packaged for shipping.

9 Disposal

The electrosurgical forceps, the packaging material and the accessories must be disposed of in accordance with the regulations and laws specific to the country in which they are used.

10 About this Instructions for Use

Throughout the period of use of the electrosurgical forceps, the Instructions for Use must be kept freely accessible for every user.

©2016 DEWIMED Medizintechnik GmbH. All rights reserved. All company names, trademarks, service marks, trade names, product names and logos appearing on this document are the property of their respective owners and should be treated accordingly.