Bipolar and Monopolar Forceps

REF 90.050.21, 90.050.25, 90.050.30
90.300.07 bis 90.715.24

IFU-19-BB-GER-ENG-IT-FR-ES Datum (date): 18.04.2018

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English Bipolar and monopolar forceps

Attention
Please read all information contained in this insert. Incorrect handling and care as well as misuse can lead to premature wear of surgical instruments or risks to patients and users.

Intended Use
Bipolar (REF 90.300.07 – 90.715.24) and Monopolar (REF 90.050.21 – 90.050.30) forceps are designed to grasp, dissect and coagulate selected tissue. They must be connected to the monopolar or bipolar output of an electrosurgical generator using a suitable cable. Only the defined parameters have to be used.

The rated supply voltage of the product is:
Bipolar forceps: 500 V_p or 300 V_p (depending on the type, see catalogue „HF-Surgery“)
Monopolar forceps: 2000 V_p
The rated accessory voltage has to be greater or equal to the maximum out peak voltage, so that the product in combination with a corresponding HF-device pursues an operating mode / attitude (see IEC/DIN EN 60601-2-2).

Appropriate connecting bipolar cables:
DEWIMED bipolar cables / flag plug, as well as DEWIMED bipolar cables / 2-pin plug; ref.no 90.190.30 – 90.208.50
DEWIMED monopolar cables: ref.no 90.220.30 – 90.221.50 / 90.224.30 – 90.225.50.

Attention: Instruments for electrosurgery should be used only by persons who have been specially trained in the use of such instruments.

Contraindications
Do not use the instrument if, in the opinion of the attending physician, the risks to the patient outweigh the benefits.

Incidents which have been reported in connection with the use of electrosurgical systems:
Unintended activation with resulting tissue injury on the wrong spot and/ or damage to the equipment.
Fire in connection with surgical drapes and other inflammable materials.
Alternating current paths leading to burns on spots where the patient or user comes into contact with components without insulation.
Explosions caused by sparks in the proximity of inflammable gases.
Perforation of organs. Sudden severe bleedings.

Use and safety instructions
The non-observance of the present use and safety instructions may lead to injuries, malfunctions or other unexpected incidents.
- When using electrosurgery in patients with pacemakers or other active implants, special requirements apply (e.g. low HF-current, patient monitoring). In any case, a cardiologist or appropriate medical specialist must be consulted.
- Before initial use and any other use, all instruments have to be completely cleaned, disinfected and sterilised and their functioning has to be checked.
- It is very important to check each surgical instrument for visible damage and wear, such as cracks, breaks or insulation defects before each use. In particular areas, such as blades, tips, notches, locking and blocking devices as well as all mobile parts, insulations and ceramic elements have to be checked carefully.
- Never use damaged instruments.
- Never use the instruments in the presence of flammable or explosive substances.
- When temporarily not in use, the instrument must be placed electrically insulated from the patient.
- Activate electrosurgical current only if the contact surfaces are visible and ensure a good contact to the tissue selected. Do not touch any other metallic instruments, trocar sleeves, optics or things like that during use.

Bipolar and Monopolar Forceps

REF 90.050.21, 90.050.25, 90.050.30
90.300.07 bis 90.715.24

IFU-19-BB-GER-ENG-IT-FR-ES Datum (date): 18.04.2018

- Observe the use and safety instructions of the manufacturer of the high-frequency surgical device.

Applies for monopolar mode of operation: Ensure correct application of the neutral electrode on the patient; otherwise, there is a danger of burns.

Reprocessing
Due to the product design, the raw materials used and the intended purpose it is not possible to determine a precise limit with regard to the maximum possible number of reprocessing cycles. The serviceable life of the instruments is determined by their function as well as by a careful handling. Instruments for electrosurgery are by nature subject to increased wear depending on the type and time of use.

Preparation and transport
Immediately after each use, clean the instruments with a soft brush under cold tap water until all visible contamination is removed. Do not use fixation agents or hot water (>40°C). Storage and transport of the instruments to the reprocessing location must be ensured in a sealed container.

Machine reprocessing

- Cleaning**
Place the instruments in a basket on the insert module or on the inserts of the MIS module and start the cleaning process.
1. Pre-rinse for 1 min. with cold water
 2. Discharging
 3. Pre-rinse for 3 min. with cold water
 4. Discharging
 5. Wash for 5 min. at 55°C with a 0.5% alkaline or at 45°C with an enzymatic cleaning agent.
 6. Discharging
 7. Neutralise for 3 min. with warm tap water (>40°C) and a neutralising agent.
 8. Discharging
 9. Rinse for 2 min. with warm tap water (>40°C).
 10. Discharging

Disinfection

Machine operated thermal disinfection has to be carried out in consideration of the national requirements with regard to the A0 value (see ISO 15883).

Drying

Dry the outside of the instruments by carrying out a drying cycle of the cleaning/disinfection machine.
If necessary, manual drying may additionally be carried out using a lintfree cloth. Dry cavities by blowing with sterile compressed air.

Manual reprocessing

Ultrasonic Pre-Cleaning

1. The instrument must be inserted in an ultrasonic bath with 0.5% enzymatic cleaning detergent. Ultrasound must be applied for 15 minutes at 40°C/104°F.
2. Remove the instrument and rinse completely with cold water to remove the cleaning detergent.

Cleaning

- Fix cleaning bath with 16 ml (1.6%) per 1 Liter of Water
1. Rinse products with cold tap water (<40°C) until all visible accumulations of dirt have been removed. Remove stuck dirt by using a soft brush.
 2. Place products in the prepared cleaning bath so that they are completely submerged: 3min
 3. Clean the instrument in the bath manually using a soft brush. Brush all surfaces several times.
 4. Rinse the products thoroughly with DI water to remove the cleaning agents without residue.

Disinfection

Place the instruments in the disinfectant bath for 12 Minutes at 20°C (Cidex OPA Solution is ready for use). Rinse the products thoroughly with fully demineralised water to remove the disinfectant without residue.

Drying

Manual drying is carried out using a lintfree cloth and, in particular, for drying cavities and channels, sterile compressed air.

Functional test and packaging

Perform visual inspection for cleanliness and integrity.
If necessary, repeat the reprocessing process until the instrument is optically clean.
Packaging has to comply with ISO 11607 and EN 868 standards for packaging for sterilised instruments.

Non-stick bipolar forceps

The polished precious metal forceps tips of the bipolar forceps may tarnish similar to silver. This does not impair function.

Bipolar forceps with irrigation

The enclosed wire insert should be always inserted in the irrigation channel, except during use and cleaning, in order to prevent clogging. The irrigation channel must be rinsed very thoroughly during cleaning. The passage has to be checked after cleaning.

Sterilisation

Sterilisation of the products with fractional pre-vacuum procedure (in accordance with ISO 13060 / ISO 17665) in consideration of the respective national requirements.
- 3 pre-vacuum phases with a pressure of at least 60 mbar.
- Heating up to a sterilisation temperature of min. 132°C and max. 137°C
- Shortest exposure time: 3 min.; at most 18 min.
- Drying time: at least 10 min.

If contamination with prions (CJK) is suspected, differing national guidelines are to be followed and longer holding times (i.e. 15 min.) may apply.

Storage

Sterilised instruments have to be stored in a dry, clean and dust-free area. The applicable national guidelines must be followed.

Repairs

Never attempt to perform repairs yourself. Service and repair work may only be performed by persons qualified and trained accordingly. For any question on these matters, please contact either the manufacturer or your medico-technical department.

Attention: Defect products must pass the complete reprocessing process before being returned for repair.

Information on the validation of the reconditioning

The following testing instructions, materials and equipment have been used for validation:

Cleaning agents (for machine use):
Neodisher FA by Dr. Weigert (alkaline)
Endozime by Ruhof (enzymatic)

Cleaning agents (manual cleaning):
Cidezyme, Enzo Enzym, detergent by Johnson&Johnson

Disinfectants (manual disinfection):
Cidex OPA, Johnson&Johnson

Neutralising agent:
Neodisher Z by Dr. Weigert

Cleaning and disinfection device:
Miele G 77365 CD
Miele insert module E 327-06
Miele MIS module E 450

For details, see report.

SMP GmbH # 01707011901 (machine cleaning)
MDS GmbH # 135196-10 (man. cleaning/disinfection)
Nelson Labs # 200432706-02 (sterilisation)
MDS GmbH Testbericht 084183-10 (sterilisation)

If the chemicals and machines described above are not available, the user has to validate the used process accordingly.

Handling

During transport, cleaning, care, sterilisation and storage, all surgical instruments should be handled with maximum care. This applies particularly to blades, fine tips and other sensitive areas.

Disposal

Disposal must be carried out in accordance with the respective applicable local and national laws and regulations.

Warranty

DEWIMED Medizintechnik GmbH exclusively supplies tested and faultless products to the customers.
All products are designed and manufactured to comply with maximum quality requirements. We refuse any liability for products which, compared to the original product, have been modified, misused or handled or used in an inexperienced way.

Explanation of symbols

	LOT-Number
	Caution: Non-sterile product
	Reference number
	Attention !
	Refer to instructions for use
	CE-Mark and registration number of the Notified Body med medical device certification, Germany
	Manufacturer Production date

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Copyright by DEWIMED Medizintechnik GmbH
SMP GmbH # 01707011901 (limpieza mecánica)
MDS GmbH # 135196-10 (limpieza/desinfección manual)
MDS GmbH Testbericht 084183-10 (esterilización)