DEUTSCH

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Copyright by DEWIMED Medizintechnik GmbH
 All rights reserved
 Subject to technical changes without notice

IFU-13-BB-GER-ENG

Datum: 21.02.2020

ENGLISH

Caution
 Please read all information contained in this insert. Incorrect handling and care as well as misuse can lead to premature wear of surgical instruments.

1 Scope

These operating instructions are applicable for DEWIMED sterile YASARGIL aneurysm clips.

| | | |
|-----------|---|------------|
| 40.560.33 | - | 40.770.13T |
| 40.780.60 | - | 40.780.60T |
| 40.782.68 | - | 40.782.68T |
| 40.783.68 | - | 40.802.13T |
| 40.810.65 | - | 40.810.65T |
| 40.820.45 | - | 40.820.45T |
| 40.830.57 | - | 40.830.57T |
| 40.900.58 | - | 40.964.55T |

2 Purpose/Function

Intended purpose

Permanent YASARGIL aneurysm clips are intended for permanent occlusion of cerebral aneurysms. Whereas temporary YASARGIL aneurysm and vessel clips are intended for temporary occlusion of blood vessels and cerebral aneurysms.

VARIANTE BOOSTER CLIPS

The YASARGIL Booster clips 40.965.00 and 40.965.00T of DEWIMED Medizintechnik GmbH may only be used to intensify the closing force of a standard YASARGIL Aneurysm Clip manufactured by DEWIMED and must not be used directly to occlude an aneurysm under any circumstances. The Titanium Booster clip must be used to intensify the closing force of YASARGIL Titanium Standard aneurysm clips. The Phynox Booster clip must be used to intensify the closing force of YASARGIL Phynox standard aneurysm clips. In order to differentiate between the materials the Titanium Booster clips are blue, the Phynox Booster clips silver in color. The notches on the Booster Clip are for correct placement on the clip that is already placed and to prevent slipping. The correct positioning of the clip must be ensured.

3 Contraindication

The permanent aneurysm clips may only be used for the permanent occlusion of cerebral aneurysms and are contraindicated for all other applications.

The temporary aneurysm clips are contraindicated for all applications except for the temporary occlusion of cerebral aneurysms or blood vessels. Temporary aneurysm clips must not be used for permanent occlusion of implantation.

4 Usage/Application

Both the permanent and temporary YASARGIL Aneurysm Clips are only intended for use by trained surgeons who are familiar with the required surgical technique and the use of these medical devices.

MINI or STANDARD clips must only be operated with the MINI or STANDARD labelled application instruments. The titanium clips and application instruments can also be identified through their different colour coding: MINI - RED & STANDARD - BLUE. Phynox aneurysm clips should only be used with the appropriate phynox application forceps. Titanium aneurysm clips should only be used with the appropriate titanium application forceps.

WARNING!
 The application of the aneurysm clips with instruments of other manufacturers and/or systems is not permitted. If the assignment (Mini & Standard) to the corresponding

application instruments is not adhered to, this can lead to functional failure or overstretching of the aneurysm clip!

5 Delivery

Permanent and temporary aneurysm clips are delivered in non-sterile condition and must be reprocessed before use (see chapter "Reprocessing of aneurysm clips"). Each package contains an aneurysm clip with a product description indicating the enclosed clip's closing force, article number (REF) and serial number (SN).

6 Reprocessing

GENERAL BASICS

The aneurysm clips are delivered in non-sterile condition and must be cleaned, disinfected and sterilized before use (cleaning and disinfection after removing the transport protection packing and sterilization after packaging). Effective cleaning and disinfection is an essential requirement for effective sterilization.

As you are responsible for the sterility of the aneurysm clips during use, please ensure

- that only sufficiently device- and product-specific validated procedures are used for the cleaning/disinfection and sterilization.
- that the used devices (disinfector, sterilizer) are maintained and checked on a regular basis and that the validated parameters are complied with in every cycle.

Please also observe the legal regulations applicable in your country and the hygiene regulations of the doctor's practice/hospital. This applies in particular to the different specifications regarding effective prior inactivation.

Aneurysm clips that have come into contact with blood, tissue or bodily fluids must be disposed of. Aneurysm clips are only intended for single use! Contaminated aneurysm clips cannot be reprocessed.

CLEANING AND DISINFECTION

Pre-treatment:

No pre-treatment is required since aneurysm clips that already came in contact with a patient or were soiled must not be used again under any circumstances.

Machine cleaning/disinfection:

(disinfector/cleaning and disinfection device)

When selecting the disinfector, it must be ensured

- that the disinfector always has proven efficacy (e.g. DGHM or FDA approval or CE marking according to DIN EN ISO 15883),
- that – if possible – a tested program for thermal disinfection (A0-value > 3000 or – for older devices – at least 5 min at 90 °C) is used (if chemical disinfection is used, there is the risk of residues of the disinfectant on the aneurysm clips),
- that the used program is suited for aneurysm clips and contains sufficient rinsing cycles,
- that only sterile or low-germ (max 10 germs/ml) and low-endotoxin (max 0.25 endotoxin units/ml) water (e.g. purified water/highly purified water) is used for rinsing,
- that the air used for drying is filtered and
- that the disinfector is maintained and checked on a regular basis.

When selecting the cleaning agent system used, it must be ensured

- that it is generally suited for cleaning aneurysm clips made of metals and plastics,
- that – unless thermal disinfection is used – a suitable disinfectant with proven efficacy (e.g. VAH/DGHM or FDA approval or CE marking) is additionally used and that it is compatible with the cleaning agent used and
- that the chemicals used are compatible with the aneurysm clips (see section "Material durability").

The concentrations specified by the manufacturer of the cleaning agent/disinfectant must be complied with.

Procedure

- 1.) Put the aneurysm clips in the tray into the disinfector. Please ensure that the aneurysm clips do not touch each other.
- 2.) Start the program.
- 3.) Take the aneurysm clips out of the disinfector after the program is finished.
- 4.) Check and pack the aneurysm clips in the tray as quickly as possible after taking them out (see the sections "Checking" and "Packaging"), if necessary after additional drying in a clean place.

The general suitability of the implants for effective machine cleaning and disinfection was confirmed by an independent certified test laboratory using the disinfector G 7836 CD (thermal disinfection, Miele & Cie. GmbH & Co., Gütersloh) and the cleaning agent Neodisher mediclean (Dr. Weigert GmbH & Co. KG, Hamburg). In doing so, the procedure described above was taken into account.

CONTROL

Inspect all aneurysm clips after cleaning or cleaning/disinfection for corrosion, damaged surfaces, chipping and contamination and discard damaged aneurysm clips. Aneurysm clips that are still contaminated must be cleaned and disinfected again. Aneurysm clips with the following features must be discarded and not used again: signs of damage, incorrect jaw positioning/misalignment, twisted components, faded color coding, contamination that cannot be removed.

PACKAGING

Please pack the trays in disposable sterilisation packaging single or double packaging) and/or sterilization containers that comply with the following requirements:

- DIN EN ISO/ANSI AAMI ISO 11607
- suitable for steam sterilisation (temperature resistance up to at least 141 °C (286 °F), sufficient steam permeability)
- sufficient protection of the aneurysm clips / sterilization packaging against mechanical damage
- regular maintenance in accordance with the manufacturer's specifications (sterilization container)


STERILIZATION

Only the sterilisation procedures listed in the following may be used for sterilisation; other sterilisation procedures are not permitted.

STEAM STERILIZATION

- fractionated vacuum procedure (with sufficient drying of the product)
- Steam steriliser according to DIN EN 13060 / DIN EN 285
- Validated according to DIN EN ISO 17665 (valid IQ/OQ (consignment) and product-specific performance evaluation)
- Maximum sterilisation temperature 138°C (280°F; plus tolerance according to DIN EN ISO 17665)
- Sterilisation time (exposure time at sterilisation temperature) at least 20 min at 121°C (250°F) or at least 3 min at 132°C (270°F)/134°C (273°F)

The general suitability of the implants for effective steam sterilization was confirmed by an independent certified test laboratory using the vapor sterilizer Systec V-150 (Systec GmbH Labor-Systemtechnik, Wetztenberg) and the fractionated vacuum procedure. In doing so, typical conditions in hospitals and doctor's practices as well as the procedure described above were taken into account.

Flash sterilization is never permissible

Do not use hot air sterilization, radio sterilization, formaldehyde or ethylene oxide sterilization or plasma sterilization, either.

STORAGE

After sterilization, the aneurysm clips must be stored dry and free from dust in the sterilization packaging.

MATERIAL DURABILITY

When selecting the cleaning agents and disinfectants, please ensure that they do not contain the following components:

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Copyright by DEWIMED Medizintechnik GmbH
 All rights reserved
 Subject to technical changes without notice

IFU-13-BB-GER-ENG

Datum: 21.02.2020

ENGLISH

- organic, mineral and oxidizing acids (minimum permissible pH value 5.5)
- strong lye (maximum permissible pH value 10.9, neutral/enzymatic or slightly alkaline cleaning agent recommended)
- organic solvents (e.g. alcohols, ether, ketones, benzine)
- oxidants (e.g. hydrogen peroxides)
- halogens (chlorine, iodine, bromine)
- aromatic/halogenated hydrocarbons

Never clean any aneurysm clips and trays with metal brushes or steel wool

All aneurysm clips and trays may only be exposed to temperatures of at most 141 °C (286 °F)!

REUSABILITY

The aneurysm clips may only be brought into contact with a patient once. Contaminated aneurysm clips may not be reprocessed and must be disposed of!

Based on the product design and the material used it is not possible to determine a precise limit with regard to the maximum number of reprocessing cycles. The service life of the medical devices is determined by their function as well as by careful handling. The reprocessing procedure of pure titanium aneurysm clips can have an impact on the visual features (blue, red and gold colour coding), i.e. the colours fade. If no clear identification of the colour coding can be given, the clips must be disposed of. Damaged and/or contaminated aneurysm clips must not be used under any circumstances.

NOTE

For patients with Creutzfeld-Jakob disease (CJD), suspected CJD or possible variations, the applicable national regulations regarding reprocessing must be complied with.

We assume no liability if these instructions are not followed.

SAFEKEEPING

To protect these high-precision instruments and delicate aneurysm clips, we recommend storing the clips in a suitable container from DEWIMED.

7 CT- & MRI-Safety

Non-clinical testing demonstrated that each YASARGIL aneurysm clip is MR conditional. A patient with the aneurysm clip can be scanned safely immediately after placement under the following conditions:

STATIC MAGNETIC FIELD

- Static magnetic field of 3-Tesla or less
- Maximum spatial gradient magnetic field of 720-Gauss/cm or less

MRI-RELATED HEATING

In non-clinical testing, these aneurysm clips produced low temperature rises during MRI performed for 15 min (i.e. per pulse sequence) in the 3 Tesla MR system (3 Tesla/128 MHz, Excite, HDx, Software 14X.M5, General Electric Healthcare, Milwaukee, WI).

The experiments showed for MRI-induced heating of these aneurysm clips at a 3 Tesla transmit/receive RF body coil MRI scanner with an average whole body SAR of 2.9 W/kg (i.e. associated with a calorimetrically measured whole body average of 2.7 W/kg) that the greatest amount of heating that occurred in association with these specific conditions was equal to or less than +1.8°C.

ARTIFACT INFORMATION

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the aneurysm clip shown below. Therefore, optimisation of MR imaging parameters may be necessary for each aneurysm clip.

8 Visual Features

PERMANENT aneurysm clips made of TITANIUM can be identified by spectral color (red or blue), TEMPORARY aneurysm clips made of TITANIUM are gold colored.

PERMANENT Aneurysm Clips made of PHYNOX are of natural colour. TEMPORARY aneurysm clips made of PHYNOX have gold plated springs and spring legs.

The exact material designation of the aneurysm clip is indicated on the packaging label.

9 Handling

Familiarize yourself with correct handling of the instruments. If necessary, let your sales representative demonstrate correct usage.

Select the correct size of aneurysm clip. Ensure that the aneurysm clip can fully encircle the aneurysm. Particular attention must be paid in the case of longer standard clip models in order to counteract a shear effect. The inner surface of the jaw of the aneurysm clip should be "filled" as completely as possible.

Bring the aneurysm clip into the sterile environment while maintaining aseptic conditions.

Examine each aneurysm clip carefully, discarding damaged or misaligned aneurysm clips. Before each implantation, check the free movement of the aneurysm clip when the jaws of the application instrument are open. If this is not present, the system must not be used and must be examined.

To avoid damages, always handle the aneurysm clips with appropriate care. Never open an aneurysm clip with your fingers and/or avoid mechanical manipulation of the aneurysm clip.

Make sure that the clip is correctly inserted into the guide grooves of the jaw of the applying forceps (see Figure 1). Do not use clip application instruments from other manufacturers.

Ease the grasps to open the jaw of the applying forceps completely. In the case of forceps with a ratchet, the ratchet must not be engaged. Make sure that the aneurysm clip lies perfectly in the guide grooves. Actuate the forceps with lock only until the lock engages (jaw closes). For complete spreading of the aneurysm clip, actuate or compress the grips on the forceps as far as they will go. The lock mechanism, if present, is released automatically. When the aneurysm clip is closed, the applying forceps (in the direction of the guide groove) can be retracted from the aneurysm clip.

If the aneurysm clip has not been correctly inserted into the application instrument or is not correctly seated in the gripper jaws of the applying forceps (see Figures 2 & 3) it may slip or be damaged if the gripper jaws close around it. Furthermore, the closing force indicated on the label can be reduced. **An incorrectly inserted aneurysm clip may be ejected from the applying forceps and may present a surgical risk.**

The correct placement of the aneurysm clip in the target tissue is especially important.

It should be ensured that the aneurysm clip was placed correctly on the neck of the aneurysm and complete occlusion of the aneurysm is guaranteed. This must be checked by the surgeon (using an appropriate method) both during and immediately after insertion.

Under certain circumstances, longer Yasargil standard aneurysm clips, straight models in particular, tend to interlock the jaws (shear effect).

The following circumstances can contribute to this:

- Application of the clip at the tip of the jaw (see Fig. 4)
- Oblique application of the clip to an aneurysm
- Application of the clip to a large aneurysm and/or severely calcified aneurysm

The article number (REF) and the serial number (SN) etched on the clip should be recorded in the patient's surgical and hospital documents.

Attention: After opening and closing the YASARGIL Aneurysm Clips more than 10 times, the nominal closing

force indicated on the packaging label may be reduced by more than 5%.

10 Possible Risks

The following serious risks were reported in connection with cerebral aneurysm clips:

Slipping of the clip, shearing of the clip, ejection of the clip from the application forceps, breakage of the clip, tearing of the aneurysm, injury to the cerebral vessels.

Apart from the risks already mentioned, additional complications can occur due to the neurosurgical procedure which cannot be directly attributed to the implant.

- Haemorrhage, secondary haemorrhage
- Wound healing abnormalities, infections
- Cerebral swelling
- Cerebral infarction
- Paralysis/sensory disturbances, functional disorders
- Epileptic seizures
- Acute hydrocephalus
- Vasospasms
- Electrolyte imbalances
- Cranial nerve injuries

Each patient must be accordingly informed by the user of the possible surgical risks when cerebral aneurysm clips are used.

DEWIMED Medizintechnik GmbH does not provide any warranty for YASARGIL Aneurysm Clips which are handled contrary to the procedure recommended here.

11 About these Instruction of Use

Throughout the period of use of surgical instruments, the Instructions for Use must be kept freely accessible for every user.

- Consult instructions for use
- Note accompanying documents
Warning: Failure to comply could result in death or injury
- Do not re-use
- Reference number
- Serial number
- Delivered in non-sterile condition
- CE-sign and shortcut of our Notified Body, mdc medical device certification GmbH, Stuttgart, Germany
- MR conditional
- Manufacturer
- Date of production
- Store in dry condition
- U.S. federal law restricts this device to be sold by or on the order of a physician only.

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Copyright by DEWIMED Medizintechnik GmbH
All rights reserved
Subject to technical changes without notice

IFU-13-BB-GER-ENG

Datum: 21.02.2020

ENGLISH