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# Instructions for use: Rotating and detachable instruments with HF-connector

## English

**Caution**  
 Please read all information contained in this insert. Incorrect handling and care as well as misuse can lead to premature wear or risks for the patient and user.

### 1 Scope

These operating instructions are applicable for rotating and detachable instruments with HF-connection (product group 2BG6) of the risk class IIb:

- REF** 26-xxx-000, 28-xxx-000, 28-000-001, 28-000-003, 28-xxx-101, 28-xxx-103, 26-xxx-101, 26-xxx-103, 28-000-100, 26-000-100, 28-000-200, 28-xxx-202, 28-xxx-204, 36-xxx-000, 37-xxx-000, 36-xxx-101, 37-xxx-101, 36-xxx-103, 37-xxx-103, 36-000-100, 37-000-100, 36-xxx-201, 36-xxx-203, 36-000-200, 36-000-001, 36-000-003, 38-xxx-100, 38-xxx-000, 38-xxx-101, 38-xxx-103, 38-xxx-201.

### 2 Maximum output voltage of the generator

Article-No.:	Umax/ Pmax
26-xxx-000, 26-xxx-101, 26-xxx-103, 26-000-100, 28-xxx-000, 28-xxx-101, 28-xxx-103, 28-xxx-202, 28-xxx-204, 28-000-001, 28-000-003, 28-000-100, 28-000-200,	1,2 kVp / 120 W
Article-No.:	Umax
36-xxx-000, 36-xxx-101, 36-xxx-103, 36-xxx-201, 36-xxx-203, 36-000-001, 36-000-003, 36-000-100, 36-000-200, 37-xxx-000, 37-xxx-101, 37-xxx-103, 37-000-100, 38-xxx-000, 38-xxx-100, 38-xxx-101, 38-xxx-103, 38-xxx-201.	2kVp

The rated accessory voltage must be greater than or equal to the maximum peak output voltage at which the product is operated in combination with an appropriate RF device, operating mode/setting (see IEC/DIN EN 60601-2-2). The insulation of the product is tested by the manufacturer for 25 reprocessing cycles. Durability in clinical practice depends on individual intraoperative use and hospital-specific reprocessing conditions.

### 3 Intended use

Rotating and detachable instruments with HF-Connection from REMA Medizintechnik GmbH are designed for coagulating of biological tissue as well as for grasping, dissecting, manipulating, cutting through and cutting of tissue, foreign objects or consumable material at endoscopic interventions.

### 4 Indication

Rotating and detachable unipolar instruments from REMA Medizintechnik GmbH are indicated at endoscopic interventions in minimal invasive surgery especially in laparoscopy and thoracoscopy.

### 5 Contraindication

The use of the rotatable and detachable unipolar instruments is contraindicated if a surgical therapy isn't indicated generally. Product specific contraindications don't exist up to now. Further contraindications:

- Local or generalized infections
- Immune suppressive therapy
- Blood clotting defect / disposition of embolism or thrombosis
- Serious systemic disease that make a survival after the surgery unlikely.
- The instruments are not allowed to be used by direct contact with heart, central circulatory system or central nervous system.

### 6 Normal use

Medical use only by medical professional staff. The instruments have to be inserted through a trocar sleeve. Rotating and detachable instruments with HF connector are connected via an RF cable (instrument side of the cable - 4mm socket) with the unipolar output of a medically approved electro-surgical generator (device side of the cable). REMA's rotating and detachable instruments with RF connection can be used with HF devices of following companies: Erbe, Martin, Berchtold, Codmann, Valleylab, Dewimed or other under comparable devices. The activation takes place by means of a foot switch.

### 7 Safty notice

#### WARNINGS AND PRECAUTIONS!

See section 2 of these Instructions for use for the maximum rated voltage of the rotating and detachable instruments with HF connector. If anything is unclear, contact the manufacturer. Before initial use, new rotating and detachable instruments with HF connector must be cleaned, disinfected and also sterilized according to a validated procedure (in this regard see section 9 "Reprocessing"). A visual inspection must be performed before each use (see section 10 "Visual Inspection"). It must be ensured that instrument-insert is inserted firmly in the handle. This must be done carefully, in order to avoid damage to the instruments and/ or injuries to the patient or surgical personnel. Rotating and detachable instruments with HF connector may be damaged if excessive force is applied.

- Danger of burning because of HF-current:**
- Use suitable neutral electrode. See instructions for use of the HF-generator.
  - Ensure correct application of the neutral electrode on the patient. Patient must lie on the whole surface of the neutral electrode.
  - Use instrument only if the isolation is undamaged.
  - Coagulation should only be performed if the contact surfaces are visible.


- It is not permissible to activate the instrument as long as it is in contact with metal objects and/or optics.
- Throughout the complete procedure, care must be taken that no flammable substances are present in the immediate vicinity, since otherwise a danger of explosion exists.
- Never use the instruments in the presence of explosive substances.
- Touch only the isolation with fingers and not the HF connector or the jaw parts.
- After the electro-surgical current has been turned off, the instrument tip may still be hot enough to cause burns.
- The high-frequency current used in electro-surgery may interfere with cardiac pacemakers and implanted heart defibrillators, and so affected patients must consult a cardiologist prior to the operation.

### 8 Assembly and Operation

Once correctly assembled and connected to power with a suitable unipolar cable, the device may be used in either the right or the left hand. To close jaws: compress (grip) handle. To open jaws: release (grip) handle. Coagulation current is activated by a footswitch that is part of your electro-surgical unit.


#### Disassembling of "Easy Change" model: 36-XXX-XXX, 37-XXX-XXX, 38-XXX-XXX:

1. Hold the instrument as shown on the picture below.
1. Press with the thumb on the pushbutton. Tube is now unlocked.
2. Pull the tube out of the handle

#### II. Hold the instrument as shown on the picture below.

1. Turn the tube left to unlock the closing-mechanism.
2. Pull the instrument-insert out of the tube.

#### Assembling of "Easy Change" model: 36-XXX-XXX, 37-XXX-XXX, 38-XXX-XXX:


- Assemble the instrument in opposite sequence.
1. By assembling the instrument please pay attention to, that the handle is opened completely. Therefore hold the instrument as shown on the picture below.
  2. Introduce the tube to the handle with a small turn until it engages.

#### Montage von "Standard Modell": 26-XXX-XXX, 28-XXX-XXX

	1. By assembling the instrument please pay attention to, that the handle is opened completely. For this purpose hold the handle as shown on the picture beside.
	2. The locking ball of the insert is mounted in the opening of the locking of the handle at an oblique angle.
	3. Close the handle. Tube is thereby drawn close to the handle.

## Instructions for use: Rotating and detachable instruments with HF-connector


4. Unscrew the connecting nut.

Disassembling of the standard model: 26-XXX-XXX, 28-XXX-XXX  
 Disassemble the instrument in opposite sequence.

### 9 Processing

Proper cleaning, inspection and maintenance will help ensure correct function of the surgical instrument. Clean, inspect and test each instrument carefully. Clean, disinfect and sterilize all instruments before surgery.

The nonsterile delivered instruments are to be taken from their packaging (for the first use). They must be placed into for cleaning / sterilization designated container / device. The detachable instruments must be cleaned in the disassembled state.

**WARNING!**

Special attentions throughout the cleaning and sterilization process require the following problem areas in the product:

- moving parts on the handle
- slots
- stops
- ends
- pipe properties

These sites are particularly difficult to access. For this reason, they must be treated with special care.

#### REPROCESSING INSTRUCTIONS:

##### Preparation at the Point of Use:

Remove gross soiling by submerge the instrument into cold water (<40°C) immediately after use. Don't use a fixating detergent or hot water (>40°C) as this can cause the fixation of residua which may influence the result of the reprocessing process.

##### Transportation:

Safe storage and transportation in a closed container to the reprocessing area to avoid any damage and contamination to the environment.

##### Preparation for Decontamination:

The devices must be reprocessed in an disassembled state.

#### AUTOMATED CLEANING PROCESS:

##### Pre-Cleaning:

**Brushing:**  
 Brush the disassembled instrument under cold water with a soft brush until all visible contamination is removed.

**Rinsing:**  
 Rinse tubes with Luer-Lock rinse port or adapter for cleaning with Luer Lock connector with the water jet pistol (static pressure above 4 bar) for a minimum time of 10 seconds. Use the water jet pistol to rinse holes, hinges, gaps and cavities of the handles and inserts as described before.

##### Automated Cleaning:

Alkaline process: Program No. 105  
 Washer/ Disinfector: Miele G 7735 CD (Miele)  
 Key Hole Surgery Rack E 450/1 (Miele)  
 Detergent: Neodisher FA, Dr. Weigert  
 Validation: SMP-Report-No.: 10109011407  
 Put the Instruments in a disassembled state on a special key hole surgery rack. Not suitable instruments must be opened and placed on an instrument tray and start the cycle.

Step:	Time (min)	Process step	Reagents	Temp. (°C)
1	1	Pre-cleaning	Tap water	cold
2		Drain		
3	3	Pre-cleaning	Tap water	cold
4		Drain		
5	5	Cleaning	Tap water - 0,5 % alkaline detergent (Neodisher FA, Dr. Weigert)	55
6		Drain		
7	3	Neutralization	Tap water	cold
8		Drain		
9	2	Intermediate-rinse	Tap water	cold

##### Automated Disinfection:

Automated Thermal Disinfection in washer/disinfector under consideration of national requirements in regards to A0-Value (see ISO 15883).

##### Drying:

Drying of outside of instrument through drying cycle of washer/disinfector. Insufflate cavities of instruments by using sterile compressed air.

##### Visual-/functional inspection:

Perform visual inspection for cleanliness and damages; assembly and functional inspection according chapter 10. If necessary, repeat the reprocessing process until the instrument is optically clean.

##### Maintenance:

After cleaning and before sterilization, treat all instruments with an oil which is considered as being physiologically safe (paraffin oil according to DAB 8 of Ph. Eur. of USP XX), especially their blades, ends, stops, snaps and all moveable parts. Distribute the lubricant uniformly in the joint by opening and closing the instrument several times. Use only instrument oils, which are approved for sterilization and taking into account the maximum sterilization temperature applied.

##### Packaging (in case of packaging before sterilization):

Packaging has to comply with ISO 11607 and EN 868 standards for packaging for sterilised instruments.

##### STERILIZATION:

Sterilization of instruments by applying a fractionated pre-vacuum process (according. ISO 13060 / ISO17665) under consideration of the respective country requirements.

Parameters for the pre-vacuum cycle:

3 prevacuum phases with at least 60 milli bar

Heat up to a minimum sterilization temperature of 132°C; maximum temperature 137°C

Minimum Holding time: 3 min (full cycle)

Drying time: minimum 1min

Validation: SMP-Report-No.: 28012

##### Packaging (in case of packaging after sterilization):

Packaging has to comply with standards for packaging for sterilised instruments.

##### Storage:

Sterilised instruments have to be stored in a dry, clean and dust-free area at moderate temperatures from 5°C to 40°C.

#### MANUAL REPROCESSING:

##### Pre-Cleaning:

Immerse the disassembled instrument into cold tap water for at least 10 minutes.

Activate all movable parts several times up to the stop.

Brush the instrument with a soft brush until all visible impurities are removed.

Flush the Instrument with a water jet pistol for minimum 20 seconds.

##### Ultrasonics:

Immerse the instrument completely into an ultrasonic bath with enzymatic detergent (0,8%) and treat with ultrasonic (frequency 35 kHz) for at least 10 minutes at 45°C.

Avoid sonic shadow.

Parameters:

- 45°C
- 10 min.
- 0,8 % cleaning solution
- 35 kHz

The instrument is taken out of the bath and flushed with a water jet pistol for minimum 20 seconds.

Ultrasonic bath: Bandelin Sonorex RK 1028 H

Detergent: Cidezime/Enzol of company ASP (enzymatic)

Validation: SMP-Report-No.: 15812

##### Manual Chemical Disinfection:

Immerse the instrument completely into a 4%-Mucocit-T solution (Company Merz Hygiene GmbH) at room temperature as specified by the manufacturer.

Validation: SMP-Report-No.: 26913

##### Drying:

Manual drying is carried out using sterile compressed air - in particular take care for drying cavities and channels.

(Continuation of the process:

##### Visual inspection

##### Maintenance

##### Packaging

##### STERILIZATION

##### Storage

as described in the above section "AUTOMATED CLEANING PROCESS".)

If the chemicals and machines described before are not available, the user is obliged to validate the process used.

## 10 Visual/ functional Inspection

**WARNING !**

It is very important to carefully examine each instrument before use for cleanliness, breaks, cracks and malfunctions of the handle, disassembling mechanism and jaw part. It is especially essential to check areas such as blades, points, ends, stops and snaps as well as all movable parts. Check insulation for cuts, voids, cracks, pressure marks, tears, abrasions, etc. Do not use unclean, damaged or not functional instruments.

## 11 Liability / Warranty Claim

REMA Medizintechnik GmbH as manufacturer and seller of these products assumes no liability for immediate or ensuing damages resulting from improper use or handling, especially due to failure to use as directed or improper processing and service.

If non-authorized persons undertake repairs or alterations of the products the manufacturer cannot be held liable; furthermore all right of guarantee are forfeited.

## Instructions for use: Rotating and detachable instruments with HF-connector

### 12 Repair and modification

Never attempt to perform repairs yourself. Service and repair work may only be performed by persons qualified and trained accordingly. Defective Surgical instruments may be sent to REMA Medizintechnik GmbH for repair. Unauthorized modifications and repairs are strictly prohibited and will entail invalidation of the manufacturer's warranty.

 **WARNING !**

**Attention:** Defect products must pass the complete reprocessing process and be signed as "hygienically safe" before being returned for repair.

### 13 Life Cycle

Due to the product design, the raw materials used and the intended purpose it is not possible to determine a precise limit with regard to the maximum possible number of reprocessing cycles. The serviceable life of the instruments is determined by their function as well as by a careful handling. Instruments for electrosurgery are by nature subject to increased wear depending on the type and time of use.

### 14 Packaging / storage / transportation / handling

Rotating and detachable instruments with HF connector must be stored in a clean and dry environment. They should be individually stored in a protective container with individual compartments or heat-sealed in film. The instruments must always be handled with the utmost care during transportation, cleaning, upkeep, sterilization and storage. This applies in particular for fine tips and other sensitive areas.

The nonsterile delivered products must be taken out of the original packaging for reprocessing.

### 15 Returns

 **WARNING !**

Returns will be accepted only if they have passed the complete reprocessing process and are marked appropriately. Returns must be securely packaged for shipping.

### 16 Disposal

Instruments and the packaging material must be disposed of in accordance with the regulations and laws specific to the country in which they are used.

### 17 About this Instruction for Use

Throughout the period of use of medical devices, the Instructions for use must be kept freely accessible for every user.

	Note accompanying documents Warning: Failure to comply could result in death or injury		Protect from moisture
	Follow instructions for use		Keep away from sunlight
<b>REF</b>	Article number (see label)	<b>LOT</b>	Batch code (see label)
	Non-sterile		Manufacturer
	CE-Mark of the notified body mdc medical device certification GmbH, Stuttgart		Do not dispose of residual waste

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Copyright by REMA Medizintechnik GmbH. Alle Rechte vorbehalten.
Copyright by REMA Medizintechnik GmbH. All rights reserved.