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# Instructions for Use

## Rigid Endoscopes

Copyright by REMA Medizintechnik GmbH  
 All rights reserved  
 Subject to technical changes without notice

IFU-31-BB-GER-ENG Date: July 06<sup>th</sup> 2021

**English**

**Products:**

- Arthroscopes** REF 19-17500 to 19-40082    **Cystoscopes** REF 31-00092 to 31-00097.14    **Hysteroscopes** REF 31-00067 to 31-00081.12
- Laparoscopes** REF 31-00010 to 31-00056    **Rhinoscopes/Sinusoscopes** REF 31-00111 to 31-00119.17    **Thoracoscopes** REF 31-00098 to 31-00099.06

**Symbols (Instructions for Use / Label / Product):**

- |  |  |   |                              |
|--|--|---|------------------------------|
|  | Attention !  |  | Observe Instructions for Use |
| <span style="border: 1px solid black; padding: 1px;">REF</span>                  | Reference number                                     | <span style="border: 1px solid black; padding: 1px;">LOT</span>                   | Batch number                 |
|  | CE-Mark and registration number of the Notified Body |  | Caution: Non-sterile product |
| 0483   |  |  | Keep away from sunlight      |
|  | Protect from moisture                                |   | Manufacturer                 |

**Attention! We deliver all products in clean, but non-sterile condition!** The instruments must be cleaned and sterilized according to below conditioning instruction before **initial use**, any further use and before return for repair.

**General Brief Advice:** The products must be used for their intended use in medical sectors only by correspondingly trained and qualified expert staff. The treating doctor or physician, purchaser or user is responsible for the selection of instruments for defined applications according to the operative use and needs to have had an appropriate training, specific information and sufficient experience for the handling of the instruments.

**Rema Medizintechnik GmbH** as the manufacturer and seller of the products rejects any responsibility for immediate damage or secondary damage caused by inappropriate use, handling or sterilization and maintenance. A repair of the instruments by firms or persons that are not authorized to repair by **Rema Medizintechnik GmbH** will cause an expiration of warranty. Before any application of the instruments these are to be checked for breaks, cracks, bendings, damages and operational reliability. The areas working end, jacket tube, optics and all movable parts must be especially carefully checked. Worn, corroded, deformed, porous or otherwise damaged products must absolutely be replaced and eliminated.

**Specified Normal Use:** The products enable the insight into the body at diverse spots. They hereby serve for diagnostics and visualization in the course of "Minimally Invasive Interventions" or through natural orifices. The responsibility for choice and combination and for the specified normal use is born by the medical scientist or user. This also includes appropriate training and instruction as well as sufficient information and experience in application of the products.

**Handling:** Handle the product carefully. Even a fall of low height or contact with other products may damage the endoscope. Avoid beats or shocks of any kind. Touch the endoscope at the housing (chamber) or at the eye funnel. The endoscope must not be deformed or bended. It is not intended for manipulation. Check connections for compatibility before connecting. Do not brute force.

**Storage:** The products have to be stored in a clean dry ambience and single in their packaging or in a protecting container with individual compartments. Protect areas such as working ends and oculars with corresponding tubes, protecting caps, gauze or fabric. Especially consider that no chemicals are stored in close proximity to the storage place.

**Materials:** The utilized materials are approved for medical technology and have been applied for years. We permanently execute extensive quality checks. Should there nevertheless be any complaint regarding the utilized materials, please refer to our customer service.

**We recommend the following procedures for the reprocessing of the endoscopes (alternatively RKI method).**

**Manual Cleaning (Pre-cleaning):**

The products preferably ought to be disinfected and cleaned immediately after use. Soiling should not dry out on the objects in order not to additionally aggravate the disinfection and cleaning. When doing so the following points need to be considered:

- The solutions which are used during the manual cleaning procedure must be applied according to the manufacturer's indications.
- A suitable brush must be used for cleaning of cannulas, blind holes and hollow parts in order to reach every spot. • Objectives and oculars (glass surfaces and glass fibre surfaces) must be cleaned with special care.
- Remove blood and other residuals with a soft brush and a mild neutral or basic (not in case of aluminium) cleaning agent. Never use metallic brushes or sponges for manual cleaning.
- In order to guarantee functionality of the products, make sure that all moving parts are thoroughly cleaned.
- Clean joints and connections in an open and closed position.
- Dismantle the products as completely as possible for reprocessing.
- Products must be stored upon easy to be cleaned carriers, such as e.g. sterilizing perforated trays/screen baskets.

**Recommendation: Bode manual instrumente detergent Bodedex forte (ph-neutral detergent).** Use a 1% solution. Water temperature of 20°C to 38°C.

Distilled water according to DIN/EN 285	Bodedex forte
1 Liter	0,01 Liter = 10 ml
2 Liter	0,02 Liter = 20ml
5 Liter	0,05 Liter = 20 ml

**Ultrasonic Treatment (accessories only)**

**Attention: Endoscopes and cables cannot be cleaned in an ultrasonic bath!**

For cleaning in ultrasonic bath the products must be placed upon sterilizing perforated trays/ screen baskets in a dismantled and opened condition. Since warm water without additives does not show satisfying cleaning results, a suitable cleaning agent must be added to the water. Manufacturer's indications regarding concentration must be considered.

The temperature of the cleaning solvents in the ultrasonic bath must be considered according to the manufacturer's indications of the utilized cleaning agent. Too important contamination will affect the cleaning result. Therefore the cleaning solution must be replaced in intervals according to the manufacturer's indications. Sonication periods must be chosen according to the manufacturer's indication of the cleaning agent. Basically ultrasonic cleaned instruments must undergo a rinse cycle. After the ultrasonic treatment the products ought to be checked for loosened parts (bolts...). Rinsing is to be effected by means of demineralized or distilled water in order to avoid water spots.

## Instructions for Use Rigid Endoscopes

Copyright by REMA Medizintechnik GmbH  
All rights reserved  
Subject to technical changes without notice

IFU-31-BB-GER-ENG Date: July 06<sup>th</sup> 2021

### English

#### Machine Cleaning:

- Cleaning and disinfection in cleaning automats is permitted after having executed a manual pre-cleaning.
- All handling and loading instructions of the machines' manufacturer have to be followed. Only the cleaning agents recommended for the specific applications by the manufacturer of the cleaning automat ought to be utilized.
- Load products into the appliance in dismantled, opened condition. Place products in a way that water can drain off from cannulas, blind holes and hollow parts.
- dismantle instruments as completely as possible for cleaning
- for instruments with long resp. narrow hollow spaces the procedures are only suitable if these hollow spaces can be floated by the hot disinfecting agent.
- During removing instruments cannulas and blind holes etc. must be checked for visible soiling. If necessary, repeat cycle or clean manually again.

#### Chemical Disinfection:

- The solvents utilized with chemical disinfection must be applied according to the manufacturer's indications of the utilized solvent.
- The applying thinning of the chemical agents is to be prepared with the use of pure water. Any addition of further cleaning agents is not allowed.
- With application of chemical agents the manufacturer's indications (residence time and concentration) must be exactly be considered.
- Disinfection solvents are to be utilized freshly daily. Multiple uses may cause the following problems:  
Increase of concentration by evaporation (risk of corrosion) or too big debris load (corrosion and decrease of effect).
- After disinfection thorough rinsing with clear floating water is basically absolutely necessary. Use of demineralized water will avoid water spots.
- The products must be sufficiently dried immediately after completion of the cleaning resp. rinsing cycles.

#### Vapour Sterilization:

Recommended sterilizing method: Vapour sterilization with fractional vacuum.

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Holding time: min. 5 minutes, max. 7 minutes

Fractionation/Pre-Vacuum = yes

The manufacturer's indications (instruction manual) of the sterilizing units are to be considered.

#### Inspection and Maintenance:

Let products cool down to ambient temperature.

Check components for damages and deformations in dismantled condition.

Products must be assembled for operational test.

You may return defected products to us in cleaned condition including proof of decontamination.

Distilled water (specifications):

Parameter	Unit	Specification
pH-value	none	6-7
Electrical conductivity	µs/cm	<0,1
Total drying residue	mg/l	1 <10
Chloride	mg/l	<0,01
Nitrate	mg/l	<0,002
Nitrat-Nitrogen	mg/l	<0,25
Sulfate	mg/l	<0,001
Magnesium	mg/l	<0,02
Spodium Nitrate	mg/l	<0,05
Nitrite	mg/l	<0,002
Phosphate	mg/l	<0,002
SiO <sub>2</sub>	mg/l	<0,01
Tin	mg/l	<0,01
Total hardness	mmol/l	<0,002