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Instructions for Use with reprocessing and safety instructions

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Unter Haßlen 14
78532 Tuttlingen

Tel.: +49(0) 7462 / 923 93 0

Fax: +49(0) 7462 / 98 93 33

Email: info@rema-surgery.de

Website: www.rema-surgery.de

Microfracturing Instruments for Arthroscopy

IFU-33-BB-GER-ENG

Date: 2021-07-06

Copyright by REMA Medizintechnik GmbH

All rights reserved

Subject to technical changes without notice

ENGLISH

Caution

Please read all information contained in this insert. Incorrect handling and care as well as misuse can lead to premature wear of surgical instruments.

1 Scope

These operating instructions are applicable for REMA reusable Microfracturing Instruments with article no. 15-18XXX-XX.

2 Intended use

Microfracturing instruments are used to stimulate cartilage replacement tissue (fibrous cartilage tissue) during arthroscopic surgery in patients with cartilage damage.

The doctor in charge respectively the user is responsible for the selection of the instruments for specific applications.

3 Indication

Locally limited cartilage damage (from stage IIb according to Noyes) of traumatic or degenerative genesis with well preserved adjacent cartilage on the knee, hip, shoulder and ankle joint.

4 Contraindication

- Leg axis malpositions with overloading of the affected joint Part
- Generalized arthrosis and corresponding cartilage bald patches of both joint partners
- Systematic arthritis and joint instability
- Patient age over 60 years (relative contraindication)

5 Normal use

The instruments were delivered nonsterile and need to be processed before using according to part 9 "processing". The instruments may only be used by specialists who are familiar with the necessary arthroscopic surgical technique for stimulating cartilage replacement tissue. Further information for a secure use of the instruments are described in part 7 "Safety notice".

6 Liability / Warranty Claim

REMA Medizintechnik GmbH as manufacturer and seller of these products assumes no liability for immediate or ensuing damages resulting from improper use or handling of the instruments, especially due to failure to use as directed or improper processing and care.

If non-authorized persons undertake repairs, alterations of the products, or combination of REMA products with the products of other manufacturer - the manufacturer cannot be held liable; furthermore all right of guarantee are forfeited.

7 Safety notice

WARNING !

The following safety aspects must be considered:

Due to the way they are constructed, the instruments have only limited stability. Excessive use of force (for example by turning or by levering) can lead to damages and influences the function of the product. Do not lever or hit microfracturing instruments with a hammer!

During transportation, cleaning/sterilization and storage, surgical instruments are to be always treated with the greatest care – especially fine tips and other sensitive parts.

Surgical instruments corrode when coming into contact with aggressive substances. This has negative effects on their function. Therefore, please follow the instructions of the cleaning agents.

Proper care is necessary for the save use of medical products.

Carry out a careful visual inspection for damage, in particular brittle or cracked rubber and plastic parts, cracks, corrosion, dirt, sharp edges and rough surfaces.

As a safety measure against potential injuries of patients and users, this check must always be performed after receipt of the instrumentation and before each use.

The instruments must be kept in clean and dry surrounding. Please ensure sufficient protection of sharp and pointed instruments.

The instrumentation is delivered non-sterile. Prior to the first time and subsequent times of use and before sending the products in for repair, the instrumentation has to be prepared according to a validated procedure (DIN EN ISO 17664) (In this regard see Section 9 "Processing").

A visual inspection must be performed before each use (see section 10 "visual inspection").

8 Life cycle

The expected life cycle of the product up to 5 years, depends of the proper use and the number of cycles.

Before each application of the product perform a visual inspection according section 10 "visual inspection".

9 Processing

Proper cleaning, inspection and maintenance will help ensure correct function of the surgical instrument. Clean, inspect and test each instrument carefully.

Clean, disinfect and sterilize all instruments before surgery.

The instruments are to be taken from their packaging (for the first use). They must be placed in a designated for cleaning / sterilization container / device.

REPROCESSING INSTRUCTIONS:

Preparation at the Point of Use:

Remove gross soiling by submerge the instrument into cold water (<40°C) immediately after use.

Don't use a fixating detergent or hot water (>40°C) as this can cause the fixation of residua which may influence the result of the reprocessing process.

Transportation:

Safe storage and transportation in a closed container to the reprocessing area to avoid any damage and contamination to the environment.

AUTOMATED CLEANING PROCESS:

Pre-Cleaning:

Brushing:

Brush the instrument under cold water with a soft brush until all visible contamination is removed.

Automated Cleaning:

Alkaline process: Program No. 105
Washer/ Disinfect- Miele G 7735 CD (Miele)
tor: Instrument Rack: E 327-06, E439 (Miele)
Key Hole Surgery Rack E 450/1 (Miele)
Detergent: Neodisher FA, Dr. Weigert
Validation: SMP-Report-No.: 10109011407

Place the instruments in a tray on the rack and start the cleaning process.

Step:	Time (min)	Process step	Reagents	Temp. (°C)
1	1	Pre-cleaning	Tap water	cold
2		Drain		
3	3	Pre-cleaning	Tap water	cold
4		Drain		
5	5	Cleaning	Tap water - 0,5 % alkaline detergent (Neodisher FA, Dr. Weigert)	55
6		Drain		
7	3	Neutralization	Tap water	cold
8		Drain		
9	2	Intermediate- rinse	Tap water	cold

Automated Disinfection:

Automated Thermal Disinfection in washer/disinfecter under consideration of national requirements in regards to A0-Value (see ISO 15883).

Drying:

Drying of outside of instrument through drying cycle of washer/disinfecter.

Visual inspection:

Perform visual inspection for cleanliness and visual inspection according chapter 10.

If necessary, repeat the reprocessing process until the instrument is optically clean.

Maintenance:

After cleaning and before sterilization, treat all instruments with an oil which is considered as being physiologically safe (paraffin oil according to DAB 8 of Ph. Eur. of USP XX). Use only instrument oils, which are approved for sterilization and taking into account the maximum sterilization temperature applied.

Packaging (in case of packaging before sterilization):

Packaging has to comply with ISO 11607 and EN 868 standards for packaging for sterilised instruments.

STERILIZATION:

Sterilization of instruments by applying a fractionated prevacuum process (according. ISO 13060 / ISO17665) under consideration of the respective country requirements.

Parameters for the pre-vacuum cycle:

- 3 prevacuum phases with at least 60 milli bar
- Heat up to a minimum sterilization temperature of 132°C; maximum temperature 137°C
- Minimum Holding time: 3 min (full cycle)
- Drying time: minimum 1min

Validation: SMP-Report-No.: 28012

Packaging (in case of packaging after sterilization):

Packaging has to comply with standards for packaging for sterilised instruments.

Storage:

Sterilised instruments have to be stored in a dry, clean and dust-free area at moderate temperatures from 5°C to 40°C.

MANUAL REPROCESSING:

Pre-Cleaning:

Immerse the instrument into cold tap water for at least 10 minutes.

Brush the instrument with a soft brush until all visible impurities are removed.

Flush the Instrument with a water jet pistol for minimum 20 seconds.

Ultrasonics:

Immerse the instrument completely into an ultrasonic bath with enzymatic detergent (0,8%) and treat with ultrasonic (frequency 35 kHz) for at least 10 minutes at 45°C. Avoid sonic shadow.

Parameters:

- 45°C
- 10 min.

Instructions for Use with reprocessing and safety instructions

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Unter Haßlen 14
78532 Tuttlingen

Tel.: +49(0) 7462 / 923 93 0

Fax: +49(0) 7462 / 98 93 33

Email: info@rema-surgery.de

Website: www.rema-surgery.de

Copyright by REMA Medizintechnik GmbH

All rights reserved

Subject to technical changes without notice

IFU-33-BB-GER-ENG

Date: 2021-07-06

ENGLISH

- 0,8 % cleaning solution
- 35 kHz

The instrument is taken out of the bath and flushed with a water jet pistol for minimum 20 seconds.

Then rinse the instrument with cold tap water again.

Ultrasonic bath: Bandelin Sonorex RK 1028 H

Detergent: Cidezyme/Enzol of company ASP (enzymatic)

Validation: SMP-Report-No.: 15812

Manual Chemical Disinfection:

Immerse the instrument completely into a 4%-Mucocit-T solution (Company Merz Hygiene GmbH) at room temperature as specified by the manufacturer.

Validation: SMP-Report-No.: 26913

Drying:

Manual drying is carried out using sterile compressed air.

(Continuation of the process:

- **Visual inspection**
- **Maintenance**
- **Packaging**
- **STERILIZATION**
- **Storage**

As described in the above section "AUTOMATED CLEANING PROCESS".)

If the chemicals and machines described before are not available, the user is obliged to validate the process used.

10 Visual Inspection

WARNING !

It is very important to inspect every surgical instrument for breaks, cracks, wear, contamination, corrosion or malfunction before each use. The tips in particular should be checked carefully.

Do not use damaged instruments.

11 Repair and modification

Never attempt to perform repairs yourself. Service and repair work may only be performed by persons qualified and trained accordingly. Defective Surgical instruments may be sent to REMA Medizintechnik GmbH for repair. Unauthorized modifications and repairs are strictly prohibited and will entail invalidation of the manufacturer's warranty.

WARNING !

Attention: Defect products must pass the complete reprocessing process and be signed as "hygienically safe" before being returned for repair.

12 Packaging, storage, transportation, handling

Surgical instruments must be stored in a clean and dry environment. They must be individually stored in a protective container with individual compartments or heat-sealed in film. The instruments must always be handled with the utmost care during transportation, cleaning, upkeep, sterilization and storage. This applies in particular for fine tips and other sensitive areas.

13 Returns

WARNING !

Returns will be accepted only if they have passed the complete reprocessing process and are marked appropriately. Returns must be securely packaged for shipping.

14 Disposal

Surgical instruments, the packaging material and the accessories must be disposed of in accordance with the regulations and laws specific to the country in which they are used.

15 About these Instructions of Use

Throughout the period of use of surgical instruments, the Instructions for Use must be kept freely accessible for every user.

	Note accompanying documents Warning: Failure to comply could result in death or injury
	Follow instructions for use
	Article number
	Batch code
	CE-Mark MP Class I
	Non-sterile
	Protect from moisture
	Keep away from sunlight
	Do not dispose of residual waste
	Manufacturer