



Fibre Light cable, autoclaveable at 2,3 bar (134°C)

Equipment adapter	Endoscope adapter	Ø mm	Length 1800 mm	Length 2300 mm	Length 3000 mm
Dewimed, Storz, Aesculap	Dewimed, Storz, Aesculap	3,5	HW31-00150	HW31-00151	HW31-00152
Dewimed, Storz, Aesculap	Dewimed, Storz, Aesculap	4,8	HW31-00200	HW31-00201	HW31-00202
Winter, Ibe	Dewimed, Storz, Aesculap	3,5	HW31-00160	HW31-00161	HW31-00162
Winter, Ibe	Dewimed, Storz, Aesculap	4,8	HW31-00210	HW31-00211	HW31-00212
Wolf, HSW	Wolf, HSW	3,5	HW31-00170	HW31-00171	HW31-00172
Wolf, HSW	Wolf, HSW	4,8	HW31-00220	HW31-00221	HW31-00222
ACM, British	ACM, British	3,5	HW31-00180	HW31-00181	HW31-00182
ACM, British	ACM, British	4,8	HW31-00230	HW31-00231	HW31-00232

High Power cables for light sources with high performance, autoclavable at 2,3 bar (134°C)

Equipment adapter	Endoscope adapter	Ø mm	Length 1800 mm	Length 2300 mm	Length 3000 mm
Dewimed, Storz, Aesculap	Dewimed, Storz, Aesculap	3,5	HW31-00540	HW31-00541	HW31-00542
Dewimed, Storz, Aesculap	Dewimed, Storz, Aesculap	4,8	HW31-00500	HW31-00501	HW31-00502
Winter, Ibe	Dewimed, Storz, Aesculap	3,5	HW31-00550	HW31-00551	HW31-00552
Winter, Ibe	Dewimed, Storz, Aesculap	4,8	HW31-00510	HW31-00511	HW31-00512
Wolf, HSW	Wolf, HSW	3,5	HW31-00560	HW31-00561	HW31-00562
Wolf, HSW	Wolf, HSW	4,8	HW31-00520	HW31-00521	HW31-00522
ACM, British	ACM, British	3,5	HW31-00570	HW31-00571	HW31-00572
ACM, British	ACM, British	4,8	HW31-00530	HW31-00531	HW31-00532